



CD/61304



CD/61104



CD/61204



CC/47022



CC/47060



CC/47062

Neutracon



A specialised surface active cleaning agent for safe cleaning and decontamination of sensitive materials which may be corroded, etched, discoloured or weakened by acidic or alkaline cleaning agents. Biodegradable, phosphate free, non-corrosive, non-flammable and totally rinsable. Specifically formulated to provide an effective near neutral cleaning agent with maximum cleaning and decontaminating properties.

Decon Neutracon is safe to use on non-ferrous metals, soft glass, coated lenses, polymer surfaces and other sensitive materials which must not be cleaned by acid or alkaline cleaning agents.

Cat. No.	Size, litre
CD/61302	1
CD/61304	5

Decon Decomatic



A liquid surface active cleaning agent for use in all types of laboratory, medical and industrial glassware washing machines. Biodegradable, bactericidal, phosphate free, non-foaming, non-flammable and totally rinsable. Decon Decomatic has been formulated to provide outstanding detergent and free rinsing action for improved results from every type of glassware washing machine which accepts a main wash detergent in liquid form. May also be used as a pre-soak prior to machine washing.

Cat. No.	Size, litre
CD/61104	5
CD/61106	10

Decon Acid Rinse



Acid rinse (neutraliser) and pre-wash (preliminary) cleaner. For use in all types of laboratory, medical and specialised industrial glassware washing machines used in conjunction with main wash detergents. Biodegradable, bactericidal, non-foaming, non-flammable and totally rinsable. Decon Acid Rinse improves the performance of glassware washing machines by preventing the carry-over of traces of the main wash detergent which may be alkaline or contain phosphates and by inhibiting the deposition of salts from 'hard' water supplies.

Cat. No.	Size, litre
CD/61204	5

Paper towels and cleaning cloths.

Paper towel, perforated roll, sheet size 24 x 46cm.

Cat. No.	Code	Description
CC/47022	7285	Single ply, blue, for removal of lighter soiling from surfaces, pack of 24 rolls
CC/47030	7308	Two ply, blue, strong, fast absorbing with good wet strength pack of 24 rolls
Accessories		
CC/47034	7041	Dispenser for single and double ply paper towels

Cleaning cloth, blue, centrefeed roll of 150 sheets highly absorbent, strong and durable (wet or dry), resists most solvents, reduces residue, allows fast light-duty cleaning.

Cat. No.	Code	Description
CC/47060	8380	Cleaning cloth roll, in free standing dispenser box

Cleaning cloth, re-useable, quarter folded, for heavy-duty tasks, highly absorbent, strong and durable (wet or dry), resists most solvents, reduces residue.

Cat. No.	Code	Description
CC/47062	8388	Cleaning cloth folded, in 4 water repellent wallets of 50 sheets