



BB/166

Weighing boats, disposable

Manufactured from high impact ant-static polystyrene to prevent powder adhesion.

Anti Static Diamond				
Cat. No.	Colour	Dim (mm)	Capacity (ml)	Pk Qty
BB/16642	Black	55 x 30	5	1000
BB/16644	Black	80 x 60	30	500
BB/16646	Black	120 x 100	100	250
BB/16648	White	55 x 30	5	1000
BB/16650	White	80 x 60	30	500
BB/16652	White	120 x 100	100	250
Anti Static Square				
Cat. No.	Colour	Dim (mm)	Capacity (ml)	Pk Qty
BB/16662	Black	44 x 44	7	1000
BB/16664	Black	80 x 80	100	1000
BB/16666	Black	140 x 140	250	1000
BB/16668	White	44 x 44	7	1000
BB/16670	White	80 x 80	100	1000
BB/16672	White	140 x 140	250	1000

Weighing Boats, Aluminium

Lightweight aluminium weighing boats, flat bottom, crimped wall for rigidity with finger tab to assist handling, FDA food-grade compliant.



BB/1670

Cat. No.	Dimensions, dia x height	Pk Qty
BB/16700	43 x 12mm	100
BB/16702	57 x 14mm	100
BB/16704	70 x 16mm	100

Weighing paper and Kjeldahl Weighing Boats

Kjeldahl Weighing Boats

- The fastest, safest, and most reliable way to transfer Kjeldahl samples.
- Dissolves residue-free in the digestion solution without influencing the analytical results in any way
- Made from very low nitrogen parchment paper without any glue or additives
- Transfer your samples without loss by simply dropping the entire weighing boat containing the sample into the acid solution in the Kjeldahl flask/digestion tube.

Parchment Paper

- Transparent and smooth
- Simplifies sample transfer
- Quantitative transfer from paper

Product	Grade	Thickness mm	Weight g/m ²
Weighing Boats, 0.07%	N 609	0.07	80
Parchment Paper, 0.05%	N B-2	0.04	43



BP/13032

Cat. No.	Code	Description	Pk Qty
BP/13032	10313032	609 Kjeldahl Weighing Boat	100
BP/17671	10347671	B-2 Sheets, 3 x 3 inch	500
BP/17672	10347672	B-2 Sheets, 4 x 4 inch	500
BP/17673	10347673	B-2 Sheets, 6 x 6 inch	500
BP/17670	10347670	B-2 Sheets, 12 x 12 inch	500