

Stirrer bars


Stirrer bars, P.T.F.E covered, cylindrical

Plain design with rounded ends.

Cat. No.	Length x o.d., mm	Pk Qty
SB/76702	8 x 3	5
SB/76704	12 x 4.5	5
SB/76706	25 x 6	5
SB/76708	30 x 6	5
SB/76710	45 x 8	5
SB/76712	50 x 8	5
SB/76714	60 x 9	5



SB/767




SB/762

Stirrer bars, P.T.F.E covered with pivot ring

For use in containers with bases that are slightly curved or uneven.

Cat. No.	Length x o.d., mm	Pk Qty
SB/76222	8 x 3	5
SB/76224	12 x 4.5	5
SB/76226	25 x 6	5
SB/76228	30 x 6	5
SB/76230	45 x 8	5
SB/76232	50 x 8	5
SB/76234	60 x 10	5



SB/763

Stirrer bars, P.T.F.E covered, disc

For use in tubes. The smallest size will create a vortex in a standard 10mm cuvette.

Cat. No.	Length x o.d., mm	Pk Qty
SB/76340	6 x 9	5
SB/76342	6 x 10	5
SB/76344	10 x 20	5
SB/76346	12 x 30	5



SB/764

Stirrer bars, P.T.F.E covered, oval

For use in round bottom flasks.

Cat. No.	Length x o.d., mm	Pk Qty
SB/76450	20 x 10	5
SB/76452	30 x 16	5
SB/76454	40 x 20	5
SB/76456	50 X 20	5

Stirrer bar sets

Boxed set of 12 octahedral bars containing 2 of each; 13x8, 15x8, 25x10, 38x10, 51x10 and 64x10mm

Cat. No.	Description
SB/76722	Stirrer bar, set of 12

Boxed set of 18 plain cylindrical bars containing 2 of each; 10x6, 15x6, 20x7, 25x8, 30x8, 40x8, 50x8, 60x10 and 80x10mm

Cat. No.	Description
SB/76724	Stirrer bar, set of 18



SB/76724

Magnetic stirrer bar retrievers, P.T.F.E. coated

A useful tool for removing stirrer bars from vessels

Cat. No.	Size
SB/76802	150 x 10mm
SB/76804	250 x 10mm
SB/76806	350 x 10mm



SB/7680