



BO/16020

## Pouring rings for media lab bottles, standard

**PYREX®**

- For use with 1516, 1519 and 1520 bottles
- Ensures drip-free pouring
- Clear polypropylene pouring ring
- Withstands autoclaving up to 140°C

Cat. No.	Code	For thread, mm	Colour	Pk Qty
BO/16020	4516/32D	32	Clear	10
BO/16022	4516/45D	45	Clear	10



BO/16024

## Pouring rings for media lab bottles, high temperature

- For use with 1516, 1519 and 1520 bottles
- Ideal for applications involving high temperatures
- Red Ethylene-tetrafluoroethylene pouring ring
- Withstands hot air sterilisation up to 180°C

Cat. No.	Code	For thread, mm	Colour	Pk Qty
BO/16024	4536/32D	32	Red	10
BO/16026	4536/45D	45	Red	10



BO/2230

## Spare caps and pouring rings for media lab bottles, polypropylene

**SCHOTT  
DURAN®**

Cat. No.	Code	For thread, mm	Pk Qty
<b>Caps</b>			
BO/22302	2923913	25	10
BO/22304	2923919	32	10
BO/22306	2923928	45	10
<b>Pouring rings</b>			
BO/22308	2924919	32	10
BO/22310	2924228	45	10



BO/126

## Spare phenolic caps

**PYREX®**

Spare caps and inserts for Pyrex screw cap, narrow neck, reagent bottles.

Cat. No.	Code	Description	Cap size, mm	Pk Qty
<b>Caps</b>				
BO/12602	4503/20	Screw cap with insert	20	20
BO/12604	4503/25	Screw cap with insert	25	20
BO/12606	4503/30	Screw cap with insert	30	20
BO/12608	4503/42	Screw cap with insert	42	10
BO/12610	4503/60	Screw cap with insert	60	5
<b>Inserts</b>				
BO/16602	4513/20M	Rubber/PTFE insert for cap	20	20
BO/16604	4513/25M	Rubber/PTFE insert for cap	25	20
BO/16606	4513/30M	Rubber/PTFE insert for cap	30	20
BO/16608	4513/42M	Rubber/PTFE insert for cap	42	10
BO/16610	4513/60	Rubber/PTFE insert for cap	60	5

NB: Insert is ethylene propylene with PTFE face for chemical resistance.