



Aspirators with stopcock, PE

Economical range



AS/96305

| Cat. No. | Shape | Capacity, litres |
|-----------------------|------------------------------|------------------|
| AS/96305 | Hexagonal | 5 |
| AS/96310 | Hexagonal | 10 |
| AS/96325 | Jerry can | 25 |
| Spare stopcock | | |
| AS/96301 | for 5 and 10 litre aspirator | |
| AS/96303 | for 25 litre aspirator | |



AS/16205

Aspirator bottles, polypropylene



- Squat form for easy location on shelves or bench top or in transit
- Moulded in handle for safe and secure carrying
- Clear moulded graduations allowing quick and easy volume assessment
- Leakproof heavy duty stopcock opens and closes with a quarter turn
- Autoclavable polypropylene for easy cleaning and sterilisation
- Caps have moulded drill guides for easy location of probes and feed-lines.
- Complimentary set of large bottles available (see BO/169 Series)

| Cat.No. | Ref. with stopcock | Capacity litres | Cap size mm | Height with cap & max width mm |
|-----------------------|---------------------------|--------------------|-------------|--------------------------------|
| AS/16205 | BNP05A | 5 | 45 | 270X220 |
| AS/16210 | BNP10A | 10 | 90 | 360X280 |
| Spare stopcock | | | | |
| AS/16201 | SCBNP01 | for BNP Aspirators | | |
| Spare cap | | | | |
| AS/16204 | For AS/16210 and AS/16220 | | | |

Notes
