

Microcentrifuge tubes, Safety cap

- Autoclavable with open cap at 121°C for 20 minutes
- Centrifugable up to 30000g
- RNase, DNase, DNA and pyrogen free
- High chemical resistance
- Graduated
- No risk of accidental opening during incubation
- Produced in a controlled clean room

Cat. No.	Capacity ml
CT/56150	0.5ml
CT/56152	1.5ml
CT/56154	2.0ml

Safety microcentrifuge tubes featuring a specially designed cap that securely snaps the tube shut for use where sample integrity is of the utmost importance.

- Autoclavable at 121°C, freezable to -90°C
- Closed cap water bath to 80°C
- Withstands centrifugation up to 30,000 g
- RNase, DNase and pyrogen free

Cat. No.	Description
CT/12050	Tube, microcentrifuge with safety cap, non-sterile, 5ml

Individual PCR Tubes

- 0.2ml and 0.5ml flat capped tubes are compatible with all current thermal cycler blocks with heated lids.
- Integral snap shut caps.
- Certified free from DNase, RNase and endotoxins.
- Natural colour (Also available in Red and Yellow, enquire for prices)



Cat. No.	Code	Description	Pk Qty
CT/56312	THER-02N	0.2ml, thin walled with domed cap	1000
CT/56314	THER-02NF	0.2ml, thin walled with flat cap	1000
CT/56318	THER-05NX	0.5ml, thin walled with flat cap	1000

Microtube rack, acrylic

Rack for 5ml Microcentrifuge tubes, 20 place, acrylic, houses all major brands of 5mL tubes. The tubes sit upright and are supported at two levels in the rack. The crystal clear colour of the rack allows for easy sample viewing.

Cat. No.	Description
TR/12000	Rack, Microcentrifuge tube, acrylic, 20 x 5ml

PCR Strip Tubes

- Compatible with all current 0.2ml 96 well cycler blocks with heated lids
- Compatible with all 8-Channel pipettors
- Certified free from DNase, RNase and endotoxins.

Cat. No.	Code	Description	Pk Qty
CT/56320	THER-08N	Strip of 8 x 0.2ml tubes, no caps	2 x 125
CT/56322	THER-ST8	Strip of caps, domed for 0.2ml tubes	2 x 125
CT/56330	THER-ST8F	Strip of caps, flat for 0.2ml tubes (RT-PCR Compatible)	2 x 125
CT/56324	THER-266	Strip of 8 x 0.2ml tubes with domed clapperboard style caps	2 x 125
CT/56326	THER-A8N	Strip of 8 x 0.2ml tubes with integral, individual hinged flat pierceable caps.	2 x 125

For suitable racks please see racks section on page 361.

Microcentrifuge tubes, screw cap.

Pre-assembled with natural cap, sterile, screw cap

Cat. No.	Code	Description
CT/56262	021-4058-510	0.5ml, conical bottom
CT/56268	021-4059-510	0.5ml, skirted base
CT/56264	021-4158-510	1.5ml, conical bottom
CT/56270	021-4159-510	1.5ml, skirted bottom
CT/56266	021-4208-510	2.0ml, conical bottom
CT/56272	021-4209-510	2.0ml, skirted bottom



CT/56152



CT/12050



CT/56318



TR/12000



CT/5632