

Active Salt software for 926 Industrial models only:

The active salt software adds to the versatility of the 926 Industrial Chloride analyser*. Originally designed for use within the cheese industry, it may be applied to any food product requiring salt analysis in a food processor's QA/ QC laboratory.


Active salt can collect data from a balance with RS232 output and the results from a Chloride analyser. One PC with one Active salt package can collect data from two analysers and two balances.

*Not suitable for use with the clinical version 926S

Advantages include:

- Advantages include: No transcription errors
- Instant product salt content calculations
- No need to achieve specific sample and diluent weights

Cat. No.	Description
CM/35002	Chloride analyser model 926 Industrial
CM/35004	Chloride analyser model 926S Clinical
CM/35034	Acid Salt software
CM/35006	Set of 3 anodes
CM/35008	Set of cathode and detector electrodes
CM/35010	Electrode polish for cleaning electrodes
CM/35012	Combined acid buffer, 500ml
CM/35014	Marked reaction beaker
CM/35016	Standard 200mg/l, 6 x 100ml
CM/35036	Standard 100mmol/l, 6 x 100ml



CM/35004

Chromatography columns**Chromatography column, with Rotaflo® stopcock**

Pyrex glass with sintered disc porosity 0 and Quickfit ground glass socket.



Cat. No.	Size, dia. x l, mm	Socket
CH/12012	10 x 200	14/23
CH/12022	20 x 300	19/26

With sintered disc porosity 0 and GP Rotaflo® stopcock


Cat. No.	Size, dia. x l, mm	Socket
CR20/30	20 x 300	19/26
CR20/50	20 x 500	19/26
CR30/50	26 x 500	24/29
CR60/50	40 x 500	40/38

Chromatography column, with Quickfit® socket and cone

Pyrex glass with sintered disc porosity 0. Joint size B24 socket and B19 cone



Cat. No.	Size, dia. x l, mm	Pk Qty
CH/12002	18 x 200	5
CH/12004	18 x 300	5



CH/1200

Chromatography papers**Chromatography papers****Whatman®**

Two grades of paper, 1 CHR for general purpose use and 3mm CHR for higher sample loading and blotting electrophoretic gels. Supplied in packs of 100 sheets or 100m reels.

Cat. No.	Code	Grade	Size, cm
Sheets			
CH/17002	3001-845	1 CHR	10 x 30
CH/17004	3001-861	1 CHR	20 x 20
CH/17006	3001-878	1 CHR	25 x 25
CH/17008	3001-917	1 CHR	46 x 57
CH/17010	3030-861	3mm CHR	20 x 20
CH/17012	3030-917	3mm CHR	46 x 57
CH/17014	3030-931	3mm CHR	58 x 68
Reels			
CH/17026	3001-604	1 CHR	1 cm wide
CH/17028	3001-614	1 CHR	2 cm wide
CH/17030	3001-640	1 CHR	3 cm wide
CH/17032	3001-652	1 CHR	4 cm wide
Other sizes available on request			



CH/170