



CD/1660



CD/80006



CD/96312

Containers

Containers, disposable, polystyrene, sterile



- All Sterilin polystyrene containers are CE marked to comply with the In-vitro diagnostic medical devices directive (98/79/EC)
- Aseptically manufactured from cytotoxicity tested FDA approved virgin polystyrene
- Containers are leak tested to BS EN14254 and most sizes have a choice of metal or plastic caps

Cat. No.	Code	Description	Materials base/cap	Pk Qty
7ml containers 50.2mm high, cap diameter 22.5mm				
CD/16602	129A	Bijou, no label	PS/PE	700
CD/16604	129B	Bijou, plain label	PS/PE	700
30ml containers 93.1mm high, cap diameter 30.1mm				
CD/16612	128A	Universal, no label	PS/PP	400
CD/16614	128B	Universal, printed label	PS/PP	400
CD/16616	128C	Universal, plain label	PS/PP	400
CD/16620	128SA	Universal, spoon, no label	PS/PE	400
CD/16622	128SB	Universal, spoon, printed label	PS/PE	400
CD/16624	128SC	Universal, spoon, plain label	PS/PE	400
60ml containers 60.7mm high, cap diameter 44.2mm				
CD/16630	125AM	Unlabelled, metal cap	PS/ME	300
CD/16632	125BM	Printed label, metal cap	PS/ME	300
CD/16634	125CM	Plain label, metal cap	PS/ME	300
60ml containers 60.7mm high, cap diameter 44.6mm				
CD/16636	125AP	Unlabelled, plastic cap	PS/PE	300
CD/16638	125BP	Printed label, plastic cap	PS/PE	300
CD/16640	125CP	Plain label, plastic cap	PS/PE	300
150ml containers 107.5mm high, cap diameter 54.9mm				
CD/16654	165A	No label, metal cap	PS/ME	120
CD/16656	165B	Printed label, metal cap	PS/ME	120
CD/16658	165C	Plain label, metal cap	PS/ME	120

Plastic containers with screw cap, aseptically produced.

Cat. No.	Description	Pk Qty
CD/80000	Container, universal, 30ml, polystyrene plastic cap, no label	400
CD/80002	Container, universal, 30ml, polystyrene plastic cap, plain label	400
CD/80004	Container, universal, 30ml, polystyrene plastic cap, patient label	400
CD/80006	Container, specimen, 60ml polypropylene, metal cap, no label	300
CD/80008	Container, specimen, 60ml polypropylene, metal cap, patient label	300
CD/80010	Container, specimen, 60ml polypropylene, metal cap, plain label	300
CD/80012	Container, specimen, 60ml polypropylene, plastic cap, no label	300
CD/80014	Container, specimen, 60ml polypropylene, plastic cap, patient label	300
CD/80016	Container, specimen, 60ml polypropylene, plastic cap, plain label	300
CD/80018	Container, specimen, 100ml polystyrene, metal cap, no label	160
CD/80020	Container, specimen, 100ml polystyrene, metal cap, plain label	160
CD/80022	Container, specimen, 100ml polystyrene, metal cap, patient label	160

Containers, multi-purpose with red screw cap

Can be used for urine samples, leak-proof. Available in both sterile and non-sterile versions.

Cat. No.	Capacity, ml	Sterile	Material	Graduation/ Labelling area	Pk Qty
CD/96300	30	No	Polypropylene	No	1000
CD/96302	30	Yes	Polypropylene	No	1000
CD/96304	60	No	Polypropylene	Yes	600
CD/96306	60	Yes	Polypropylene	Yes	500
CD/96308	120	No	Polypropylene	Yes	500
CD/96310	120	Yes	Polypropylene	Yes	250
CD/96312	200	No	Polypropylene	No	200
CD/96314	200	Yes	Polypropylene	No	150
CD/96316	200	Yes	Polystyrene	No	150