

Crucibles

Crucibles, Nickel, 0.5mm wall thickness

Supplied individually.

Cat. No.	Capacity	Diameter	Height
CR/94102	10ml	25mm	30mm
CR/94106	30ml	40mm	40mm
CR/94108	50ml	45mm	45mm
CR/94110	70ml	50mm	50mm

Nickel crucibles with 1.0mm, 1.5mm and 2.0mm wall thickness are available, details on request

Crucible lids for above, Nickel

Supplied individually.

Cat. No.	To fit crucible capacity
CR/94132	10ml
CR/94136	30ml
CR/94138	50ml
CR/94140	70ml

Crucibles, Zirconium, 1.0mm wall thickness

Supplied individually.

Cat. No.	Capacity	Diameter
CR/94112	10ml	25mm
CR/94114	25ml	35mm
CR/94116	30ml	40mm
CR/94118	50ml	45mm
CR/94120	70ml	50mm

Crucible lids for above, Zirconium

Supplied individually.

Cat. No.	To fit crucible capacity
CR/94142	10ml
CR/94144	25ml
CR/94146	30ml
CR/94148	50ml
CR/94150	70ml

Crucibles, stainless steel

Lower cost alternative, more durable and longer lasting than porcelain. Easier to clean and are shatter-proof. Thickness 0.7mm.

Cat. No.	Capacity, ml	Diameter,mm	Height,mm	Pk Qty
CR/94160	15	30	30	10
CR/94162	25	35	35	10
CR/94164	50	45	45	10
CR/94166	70	50	50	10

Crucible lids for above, Stainless steel

Cat. No.	To fit crucible capacity	Pk Qty
CR/94170	15ml	10
CR/94172	25ml	10
CR/94174	50ml	10
CR/94176	70ml	10



CR/9410



CR/9411



CR/9416



CR/9417