



CR/942

Crucibles, Platinum, with plain rim

All weights are approximate, tolerance of $\pm 5\%$

Cat. No.	Capacity,ml	Diameter,mm	Height,mm	Weight, gm
CR/94202	10	26	29	8.5
CR/94204	15	29	33	11.5
CR/94206	20	32.5	36	18
CR/94208	25	35	39	19
CR/94210	30	37.5	41	24
CR/94212	40	41.5	46	31
CR/94214	50	44	49	42



CR/98008

Crucible lids for above, Platinum

All weights are approximate, tolerance of $\pm 5\%$

Cat. No.	To fit crucible capacity	Weight, g
CR/94222	10ml	2.8
CR/94224	15ml	3.6
CR/94226	20ml	4.5
CR/94228	25ml	5.2
CR/94230	30ml	6
CR/94232	40ml	7.2
CR/94234	50ml	8



CR/98016

Crucible, porcelain

Glazed inside and out except base, without lids, packs of 10

Cat. No.	Capacity,ml	Diameter,mm	Height,mm
Low form			
CR/98002	5	30	19
CR/98004	10	35	22
CR/98006	17	40	25
CR/98008	21	45	28
CR/98010	34	50	32
CR/98012	63	60	38
CR/98014	91	70	44
Medium form			
CR/98016	12	35	28
CR/98018	20	40	32
CR/98020	30	45	36
CR/98022	45	50	40
CR/98024	80	60	48
CR/98026	120	70	56
Tall form			
CR/98028	15	30	38
CR/98030	26	35	44
CR/98032	35	40	50
CR/98034	50	45	56
CR/98036	72	50	62
CR/98038	130	60	75



CR/98032

Crucible lids for above, Porcelain, Pack of 10



CR/9804

Cat. No.	Diameter,mm
CR/98040	30
CR/98042	35
CR/98044	40
CR/98046	45
CR/98048	50
CR/98050	60
CR/98052	70