




CR/36512



CR/36502



CR/10400



CR/10420



CR/10450

Cryogenic storage

Cryogenic vial holder 25-place

NALGENE®

Polycarbonate. Allows single-handed vial manipulation. Bottom of well locks with base of vial. Spacing around each hold is greater than 5030-0510 to grasp vials easily. Accepts all Nalgene cryogenic vials listed. Autoclavable.

Cat. No.	Code	Capacity	Dimensions, WxDxH, mm
CR/36512	5030-0505	25 vials	97x102x22

Cryogenic vial holder 50-place

NALGENE®

White Polycarbonate. Allows single-handed vial manipulation. Bottom of each well interlocks with base of vial. Accommodates up to 50 vials in a 5 x 10 array. Accepts all Nalgene cryogenic vials listed.

Cat. No.	Code	Dimensions, WxDxH, mm
CR/36502	5030-0510	102x197x28

Storage boxes, Cryogenic, 81 and 100 place

Durable Polypropylene storage boxes are the perfect solution for freezer storage and organization. Matte surface on box sides are suitable for printing and or labelling. Boxes withstand temperatures from -80°C to +121°C and will not degrade like cardboard boxes.

81 Place boxes:

- Tight fitting lid, samples will not fall out; even when held upside down.
- Fits in standard 5¼ square, stainless steel, freezer racks.
- Holds 0.35, 1.5 and 2ml microtubes as well as most brands of vials.
- Alphanumeric grid for sample identification.
- Clear top for visual inspection, base is available in 4 fluorescent colours as well as natural or in a pack of assorted colours.
- Fits most standard 2 stainless steel freezer racks.

100 Place boxes:

- Lock-and-seal triple hinged locking front clasp keep samples secure.
- Holds 0.35, 1.5 and 2ml microtubes as well as most brands of vials.
- Dual alphanumeric grid (Inside and outside) for sample identification.
- Clear top for visual inspection, base is available in 4 fluorescent colours as well as natural or in a pack of assorted colours
- An all-black box is also available for light-sensitive samples.

Cat. No.	Description
Polypropylene 81 place Cryogenic storage box, 5 1/2 x 5 1/2 x 1 3/4 inches, pack of 5	
CR/10400	Storage box, Cryogenic 81 place, polypropylene, natural
CR/10402	Storage box, Cryogenic 81 place, polypropylene, blue
CR/10404	Storage box, Cryogenic 81 place, polypropylene, green
CR/10406	Storage box, Cryogenic 81 place, polypropylene, purple
CR/10408	Storage box, Cryogenic 81 place, polypropylene, orange
CR/10410	Storage box, Cryogenic 81 place, polypropylene, assorted colours
Polypropylene 100 place Cryogenic storage box, 5 1/2 x 5 1/2 x 2 11/32 inches, pack of 5	
CR/10420	Storage box, Cryogenic 100 place, polypropylene, natural
CR/10422	Storage box, Cryogenic 100 place, polypropylene, blue
CR/10424	Storage box, Cryogenic 100 place, polypropylene, green
CR/10426	Storage box, Cryogenic 100 place, polypropylene, purple
CR/10428	Storage box, Cryogenic 100 place, polypropylene, orange
CR/10430	Storage box, Cryogenic 100 place, polypropylene, assorted colours
CR/10432	Storage box, Cryogenic 100 place, all black for light sensitive samples

Storage box, gridless, clear

Polypropylene box without dividers, 2 inches tall, to store larger tubes (lying flat) and laboratory vials that won't fit in standard freezer boxes. The bottom of the box includes venting holes.

Cat. No.	Description
CR/10450	Storage box, gridless, clear polypropylene, 5 1/4 x 5 1/4 x 2 inches