

Storage boxes, freezer, Polycarbonate

Polycarbonate freezer boxes, with clear lid, available with 25, 81, and 100-place divisions with locator grids to facilitate sample identification. The 25 & 81 place boxes have individual cell compartment dividers and the 100-place box holds tubes in place with a pegged grid divider. Designed to withstand temperatures from -196°C to 121°C, the boxes have drain holes for venting and safety and a keyed lid prevents misalignment of the grid.

Cat. No.	Description
CR/10460	Storage boxes, freezer, 3 x 3 x 2 1/16 inches, 25 place
CR/10462	Storage boxes, freezer, 5 1/4 x 5 1/4 x 2 inches 81 place
CR/10464	Storage boxes, freezer, 5 1/4 x 5 1/4 x 3 11/16 inches 81 place
CR/10466	Storage boxes, freezer, 5 1/4 x 5 1/4 x 2 inches, 100 place



CR/10460

Cryobox

Made from strong polycarbonate, these durable cryogenic storage boxes are designed to be used at temperatures between -196°C and +121°C and are autoclavable at 120°C, 15psi (1 bar) for 20 minutes. Different styles are available to accommodate 25, 81 or 100 cryovials. A transparent cover allows the user to see the contents of the box and is keyed to the base to prevent misalignment control. Boxes are colour coded, either with coloured grids on the 25 and 81 place boxes, or coloured bases for the 100 place boxes. An innovative vial picker is supplied with each box to facilitate easy removal of vials.

25 Place Cryostore box from 1-2ml Cryovials

Cat. No.	Description	Capacity	Dimensions LxWxH,mm	Colour	Pk Qty
CR/91000	Cryostore™ storage box	25	76x76x52	Blue	8
CR/91002	Cryostore™ storage box	25	76x76x52	Green	8
CR/91004	Cryostore™ storage box	25	76x76x52	Red	8
CR/91006	Cryostore™ storage box	25	76x76x52	Yellow	8

81 Place Cryostore box for 1-2ml Cryovials

Cat. No.	Description	Capacity	Dimensions LxWxH,mm	Colour	Pk Qty
CR/91020	Cryostore™ storage box	81	133x133x52mm	Blue	4
CR/91022	Cryostore™ storage box	81	133x133x52mm	Green	4
CR/91024	Cryostore™ storage box	81	133x133x52mm	Red	4
CR/91026	Cryostore™ storage box	81	133x133x52mm	Yellow	4

100 Place Cryostore box for 1-2ml Cryovials

Cat. No.	Description	Capacity	Dimensions LxWxH,mm	Colour	Pk Qty
CR/91040	Cryostore™ storage box	100	133x133x52mm	Blue	4
CR/91042	Cryostore™ storage box	100	133x133x52mm	Green	4
CR/91044	Cryostore™ storage box	100	133x133x52mm	Red	4
CR/91046	Cryostore™ storage box	100	133x133x52mm	Yellow	4



CR/91004

Storage boxes, Cardboard, White, 133mm x 133mm

- Temperature-resistant down to -80°C
- Also available with a waterproof plastic coating
- Wide variety of inserts available (supplied separately)
- Four additional colours are available, blue, red, green and yellow
- Five different heights

Water repellent standard coating, without inserts

Cat. No.	Height	Pk Qty
CR/27010	32mm	10
CR/27012	50mm	10
CR/27014	75mm	10
CR/27016	100mm	10
CR/27018	130mm	10



CR/91020