



AC/161



AC/1630



AC/80004



AC/96102

Autoclave Accessories



Autoclave bags and stands

Strong disposable bags with blue biohazard labelling designed for contaminated waste disposal in autoclaves or incinerators. Convenient 'tissue box' cartons dispense bags individually and wire and cardboard stands are designed to hold bags open while in use.

Bags

Cat. No	Sterilin Code	Max Temp (°C)	Material	Length x Width (mm)	Pk Qty
AC/16102	509	121 (low)	HDPE	305 x 660	200
AC/16104	510	121 (low)	HDPE	305 x 660	500
AC/16106	509L	121 (low)	HDPE	406 x 610	200
AC/16108	510L	121 (low)	HDPE	406 x 610	500
AC/16110	511	121 (low)	HDPE	610 x 810	200
AC/16202	509HT	134 (high)	PP	307 x 660	200
AC/16204	510HT	134 (high)	PP	307 x 660	500
AC/16206	509LHT	134 (high)	PP	406 x 610	200
AC/16208	510LHT	134 (high)	PP	406 x 610	500
AC/16210	511HT	134 (high)	PP	610 x 810	200

Stands

Cat. No	Sterilin Code	Material	Description	Pk Qty
AC/17002	S23B	Cardboard	Holder for 511	10
AC/14002	S23C	Coated Wire	Holder for 509, 510	1
AC/14004	S23E	Coated Wire	Holder for 511	1

Bags

Biohazard autoclave bags with warning instruction and sterilization indicator field, strong 50 m Polypropylene film with broad seam at bottom. Autoclavable at 134°C (When autoclaving, the bag must not be completely sealed)

Cat No.	Dimensions, wxl (mm)	Volume (approx) Litres	Pk qty
AC/16300	300 x 500	8	500
AC/16302	400 x 780	38	500
AC/16304	600 x 800	60	500
AC/16306	700 x 1100	110	350

Autoclave Deodoriser Capsules

Anabac® fragrant capsules that release a discreet and pleasant perfumed odour, placed directly in the autoclave, they melt without residue and have no effect on the disinfection process. A variety of fragrances are available, the capsules are 100% biodegradable and are suitable for use with any autoclave, supplied in packs of 100 capsules.

Cat. No.	Description
AC/80002	Autoclave deodorant based on Eucalyptus extract
AC/80004	Autoclave deodorant based on Apple extract
AC/80006	Autoclave deodorant based on Lemon extract
AC/80008	Autoclave deodorant based on Flower extract
AC/80010	Autoclave deodorant based on Peach extract

Autoclave tape

Indicator tape for steam sterilisation.

Cat. No.	Description	Size
AC/97002	Autoclave tape	19mm x 50m

Autoclave adhesive tape (Class 4) Can be written on, for use in standard autoclaves. Irreversible colour transition from Green to Orange after 15 minutes at 121°C. Provides clean removal from autoclaved object (leaves no residue). Pack of 12 rolls.

Cat. No.	Code	Description	Size
AC/96102	TAP-SCS-E	Autoclave tape	18mm x 55m