

Cylinders, measuring, moulded graduations, PMP

- Clear, easy to read moulded graduations
- Conforms to class B tolerances
- Transparent, robust PMP
- Excellent chemical resistance
- Conforms to ISO 6706
- Large pouring spout for easy dispensing
- Easy to read as liquids have virtually no meniscus in plastic cylinders

Azlon

Cat. No.	Code	Size, ml	O.D. x Height, mm	Tolerance, ± ml	Pk Qty
CY/16828	CT0010P	10	14 x 140	0.2	10
CY/16834	CT0025P	25	20 x 170	0.5	5
CY/16840	CT0050P	50	26 x 200	1	5
CY/16846	CT0100P	100	31 x 250	1	5
CY/16856	CT0250P	250	42 x 315	2	1
CY/16866	CT0500P	500	55 x 360	5	1
CY/16876	CT1000P	1000	66 x 440	10	1
CY/16880	CT2000P	2000	83 x 535	20	

Note: Temperatures above 80°C may affect accuracy



CY/168

Cylinders, measuring, moulded graduations, PMP

Kartell

Conforms to ISO 6706 - 1981 (E) and BS 5404 Part 2 1977. Crystal clear. Autoclavable. Can be used for liquids up to 170°C. No wetting. No meniscus. Permanent graduations are moulded in for reproducibility. Chemically cleaner than glass. Chemically non-absorbent. Excellent chemical resistance. Pentagon base for added stability.

Cat. No.	Code	Capacity, ml	Grads, ml	Pk Qty
CY/46828	1570	10	0.2	10
CY/46834	1571	25	0.5	5
CY/46840	1572	50	1.0	5
CY/46846	1573	100	1.0	5
CY/46856	1574	250	2.0	1
CY/46866	1575	500	5.0	1
CY/46876	1576	1000	10.0	1
CY/46880	1577	2000	20.0	1



CY/468

Cylinders, Nessler, soda glass

For colour comparison. Supplied in matched pairs.

volac

Cat. No.	Description	Pk Qty
CN/32002	50ml single mark	pair
CN/32004	100/50ml two marks	pair



CN/32004

Cylinders, receivers, Crow pattern, Type 1, Class B

Cat. No.	Capacity, ml
CY/52202	100



CY/52202