



DE/1650




DE/36502



DE/32504



DE/3600



DE/16606

Desiccator, vacuum, polycarbonate

Azlon

- Tough transparent polycarbonate cover and base.
- Polycarbonate stopcock with PTFE plug.
- Will safely withstand high vacuum pressures making it safe and virtually unbreakable.
- Ceramic plate permits the use of incandescent crucibles.

Cat. No.	Code	I.D,mm	Plate Size, mm	Overall height,mm
DE/16502	DWA215	215	205	262
DE/16504	DWA245	245	236	296

Desiccator, vacuum, polycarbonate

NALGENE

Transparent light weight and unbreakable. For vacuum and non-vacuum applications. Holds vacuum of 711mmHg for 24 hours. For room temperature use only. Silicone O-ring seal does not require grease. Not suitable for use with organic solvents. Supplied without disc.

Cat. No.	Code	Size top I.D., mm
DE/36502	5311-0250	251

Desiccator plate, ceramic-metal

NALGENE

Unbreakable, with greater resistance to thermal shock than porcelain. Fire polished, corrosion resistant, inert non-stick glass surface bonded to metal. The plate is marked with numbered quadrants for easy location of contents. Autoclavable.

Cat. No.	Code	Diameter, mm
DE/32504	5312-0230	230

Desiccator plates, porcelain

Cat. No.	Description
DE/62030	Porcelain plate for 100mm desiccator DIN, 90 mm
DE/62032	Porcelain plate for 150mm desiccator DIN, 140 mm
DE/62034	Porcelain plate for 200mm desiccator DIN, 190 mm
DE/62036	Porcelain plate for 250mm desiccator DIN, 235 mm
DE/62038	Porcelain plate for 300mm desiccator DIN, 280 mm

Desiccator cabinets, acrylic

NALGENE

Provide dust and moisture-proof storage in refrigerators, freezers and on the bench top. Useful to store anhydrous, hygroscopic biological and chemical compounds, electronic assemblies or any components requiring controlled storage conditions. Neoprene gasket on door ensures a tight seal. Door is secured with two stainless steel latches and full height hinge. Removable lipped tray holds solid desiccant. Shelves are removable, adjustable and vented to maximise air flow.

Cat. No.	Code	Dim.(HxWxD), mm	Shelves
DE/36002	5317-0070	178x303x305	2
DE/36004	5317-0120	305x305x305	3
DE/36006	5317-0180	457x305x305	4

Desiccator cabinet, polystyrene

Azlon

- Clear, high impact polystyrene cabinet
- Space saving stacking feature
- Airtight door with foam rubber gasket
- Supplied with sealed jar of silica gel

Cat. No.	Code	Internal LxWxH, mm	External LxWxH, mm	Door opening, mm
DE/16606	DWA301	201x180x130	225x200x168	188x110

For Dewars, please see Cryogenics page 138.