



BB/167

Dishes, Aluminium

Dishes, Aluminium

Lightweight aluminium weighing boats, flat bottom, crimped wall for rigidity with finger tab to assist handling, FDA food-grade compliant.

Cat. No.	Dimensions, dia x height	Pk Qty
BB/16700	43 x 12mm	100
BB/16702	57 x 14mm	100
BB/16704	70 x 16mm	100

Dispensers

Dispensers, bottle-top, BibbyPet



The BibbyPet is a new bottle-top dispenser with an innovative design concept that makes it a real low cost alternative for routine dispensing. Based on a replaceable polypropylene/polyethylene dispensing cartridge the BibbyPet is ideal for use with aqueous solutions, alkalis, dilute acids and polar solvents such as ethanol or acetone. Both the dispensing cartridge and valves can be quickly and easily replaced. The BibbyPet has a 45mm thread and is supplied complete with adapters for bottles with 32mm and 40mm threads. Both the 10ml and the 25ml versions have an accuracy of $\pm 1.2\%$ and a repeatability of 0.2%

- Economical
- Semi-automatic
- Replaceable dispensing cartridge and valves



DI/13002

Cat. No.	Code	Description
DI/13002	BD10	10ml capacity BibbyPet with 0.2ml setting
DI/13004	BD25	25ml capacity BibbyPet with 0.5ml setting
DI/13006	BD1010	Dispensing cartridge 10ml, pack of 3
DI/13008	BD1012	Dispensing cartridge 25ml, pack of 3
DI/13010	BD1014	Discharge tube
DI/13012	BD1016	Suction tube, 250mm, pack of 2
DI/13016	BD1020	Valve set, including seals

Dispensers, bottle top, Pressmatic



The high specification and innovative design of the Pressmatics make them the ideal dispenser for almost any solution, including concentrated acids and alkali reagents. A new recirculation system cuts down on solution wastage during priming. The Pressmatic fits directly onto a 45mm thread and is supplied with three adapters to fit most laboratory bottles. Every unit is individually calibrated to 0.7% accuracy and supplied with a certificate.

- New design with recirculation valve
- Chemically resistant and non-sticking ceramic piston
- Simple & accurate volume setting
- Fully rotational on bottle for optimum positioning
- Accuracy within 0.7%
- Can be autoclaved at 121°C



DI/131 with accessory bottle

Cat. No.	Code	Capacity, ml
DI/13102	PD1R	0.2 to 1.0
DI/13106	PD5R	1.0 to 5.0
DI/13108	PD10R	2.0 to 10.0
DI/13110	PD30R	5.0 to 30.0
DI/13112	PD60R	10.0 to 60.0

Accessory bottles for Pressmatic dispensers

Cat. No.	Description
DI/13116	1000ml, amber
DI/13120	1000ml, clear