

Scalpel handles, acrylic

Swann-Morton®

Type 5A and 6A.

5A accepts standard scalpel blades style 10, 10A, 11, 12, and 15.

6A accepts standard scalpel blades style 20, 21, 22, 22A, 23, 24, 25, 26 and 27.

Cat. No.	Type	Pk Qty
DI/26002	5A	Each
DI/26004	6A	Each



DI/26002



DI/26004

Scalpel blades, carbon steel, No. 3 Fitment

Swann-Morton®

To fit handles 3, 3L, 5A and 7

Cat. No.	Code	Style	Pk Qty
Non Sterile, inner foil packs of 5			
DI/24040	0101	No.10	100
DI/24042	0102	No.10A	100
DI/24044	0103	No.11	100
DI/24046	0104	No.12	100
DI/24048	0105	No.15	100
Sterile, inner individual foil packs			
DI/24140	0201	No.10	100
DI/24142	0202	No.10A	100
DI/24144	0203	No.11	100
DI/24146	0204	No.12	100
DI/24148	0205	No.15	100



DI/24040



DI/24048



DI/24046

Scalpel blades, carbon steel, No. 4 fitment

Swann-Morton®

To fit handles 4, 4L, and 6A

Cat. No.	Code	Style	Pk Qty
Non Sterile, inner foil packs of 5			
DI/24050	0106	No. 20	100
DI/24052	0107	No. 21	100
DI/24054	0108	No. 22	100
DI/24056	0109	No. 22A	100
DI/24058	0110	No. 23	100
DI/24060	0111	No. 24	100
DI/24062	0112	No. 25	100
DI/24064	0113	No. 26	100
DI/24066	0114	No. 27	100
Sterile, inner individual foil packs			
DI/24150	0206	No. 20	100
DI/24152	0207	No. 21	100
DI/24154	0208	No. 22	100
DI/24156	0209	No. 22A	100
DI/24158	0210	No. 23	100
DI/24160	0211	No. 24	100
DI/24162	0212	No. 25	100
DI/24164	0213	No. 26	100
DI/24166	0214	No. 27	100



DI/24050



DI/24058




DI/24066

Scalpel blade remover

For the simple and safe removal of scalpel blades. Single use, plastic, sterile, disposable. For removal of all blades except size 12. Pack of 50.

Cat. No.
DI/24068



DI/24068