

Drying/warming Cabinets



Made of steel and finished in white powder coated paint with heat proof toughened glass doors and wire mesh shelves. Heating elements are waterproof and are located under the perforated floor which is finished in polished stainless steel. Temperature control is by Simmerstat giving any desired temperature up to a maximum of 50°C above ambient temperature. Smaller models have sliding doors whilst 255L models and above have hinged doors. Electrical supply, all models 220/240V 50Hz.



DC/75004

Cat. No.	Model	Fan system	Internal dims (HxWxD)	External dims (HxWxD)	Capacity	Weight	No. of shelves	Electrical loading
DC/75002	SS	None	430x750x340mm	530x760x360mm	113L	22Kg	3	500W
DC/75004	LS	None	530x990x390mm	630x1000x410mm	226L	36Kg	3	750W
DC/75022	B2	None	685x630x590mm	864x635x610mm	255L	45Kg	1	1000W
DC/75024	BC2	Circulating	685x630x590mm	964x635x610mm	255L	45Kg	1	1000W
DC/75026	BX2	Extract	685x630x590mm	1064x635x610mm	255L	45Kg	1	1000W
DC/75028	B1	None	685x1165x610mm	864x1168x610mm	480L	63Kg	1	1500W
DC/75030	BC1	Circulating	685x1165x610mm	964x1168x610mm	480L	63Kg	1	1500W
DC/75032	BX1	Extract	685x1165x610mm	1064x1168x610mm	480L	63Kg	1	1500W
DC/75006	F2	None	1430x630x590mm	1626x635x610mm	540L	75Kg	4	1500W
DC/75008	FC2	Circulating	1430x630x590mm	1726x635x610mm	540L	75Kg	4	1500W
DC/75010	FX2	Extract	1430x630x590mm	1826x635x610mm	540L	75Kg	4	1500W
DC/75014	F1	None	1430x1165x590mm	1626x1168x610mm	1000L	125Kg	4	2500W
DC/75016	FC1	Circulating	1430x1165x590mm	1726x1168x610mm	1000L	125Kg	4	2500W
DC/75018	FX1	Extract	1430x1165x590mm	1826x1168x610mm	1000L	125Kg	4	2500W
DC/75020	FCX1	Circulating/ Extract	1430x1165x590mm	1800x1168x610mm	1000L	125Kg	4	2500W



DC/75016



DC/75006



DC/75002 and DC/75004