

Bags, polythene, resealable

Polythene bags, selfseal, with write on panels (Sold in packs of 100)

Cat. No.	Approx. size (w x d), mm
BA/96430	76 x 83
BA/96432	89 x 114
BA/96433	114 x 114
BA/96434	102 x 141
BA/96436	127 x 191
BA/96437	140 x 140
BA/96438	152 x 229
BA/96439	191 x 191
BA/96440	203 x 280
BA/96442	229 x 324
BA/96444	254 x 356
BA/96446	279 x 406
BA/96448	330 x 457
BA/96449	381 x 508



BA/96430

Bags, polythene, resealable

Polythene bags, plain (Sold in packs of 1000)

Cat. No.	Approx. size (w x d), mm
BA/96470	76 x 83mm
BA/96472	89 x 114mm
BA/96474	102 x 141mm
BA/96476	127 x 191mm
BA/96478 1	152 x 227mm

Other sizes available on request



BA/96470

Bags, sample, Whirl-Pak®, sterile, write on panel

Made from polythene with tab and wire frame at the aperture allowing the bag to be opened wide. Sealing is achieved by pulling the tabs and spinning the bag around the wire. The bags can be opened and closed as required. Suitable for solids and liquids, with write on panel. Other sizes available on request.

Cat. No.	Description	Capacity ml	Pk Qty
BA/96450	Whirl pack bag, 185 x 75mm panel	118	500
BA/96454	Whirl pack bag, 230 x 115mm panel	532	500
BA/96456	Whirl pack bag, 230 x 120mm panel	710	500
BA/96458	Whirl pack bag, 380 x 190mm panel	2041	500
BA/96480	Whirl pack bag, 380 x 254mm panel	2721	250
BA/96482	Whirl pack bag, 254 x 508mm panel	3637	250
BA/96484	Whirl pack bag, 508 x 380mm panel	5441	100
BA/96486	Whirl pack bag, 230 x 115mm self-standing	532	500
BA/96488	Whirl pack bag, 380 x 125mm self-standing	1065	250
BA/96490	Whirl pack bag, 380 x 190mm self-standing	2041	250



BA/96454



BA/96480