

Multi Sub Maxi horizontal gel unit suitable for RFLP analysis, Southern and Northern blotting preps and high throughput analysis with up to 450 samples

Technical Specifications	
Active gel size	20 x 10cm + 20 x 20cm
Sample capacity	1-200 x (20cm x 10cm) 1-450 x (20cm x 20cm)
Buffer tank volume	1200ml
Combs available	
Thickness	0.75, 1.0, 1.5, 2.0mm
Number of teeth	1, 2, 4, 10, 16, 20MC, 25, 30, 36, 40MC, 50 25
Unit dimensions (WxLxH)	23 x 39.5 x 9cm

Cat. No.	Code	Description
EP/70110	MSMAXIDUO	Multi Sub Maxi, horizontal electrophoresis unit with 20 x 10 and 20 x 20cm UV tray, 2 x 20 sample combs, lid, power cables, loading guides and dams
EP/70120	MS20-LG	Loading guides
EP/70122	M20-UVS	Gel scoop, UV, 20cm
EP/70130	MS20-16-1	Comb, 16 sample, 1mm thick
EP/70132	MS20-20MC-1	Comb, 20 sample, MC, 1mm thick
EP/70134	MS20-25-1	Comb, 25 sample, 1mm thick
EP/70136	MS20-30-1	Comb, 30 sample, 1mm thick
Accessories		
EP/70140	CSL-GLT	Gel levelling table
EP/70142	CSL-RULER	Fluorescent ruler, UV



EP/70110

Vertical electrophoresis gel units

OmniPage Mini vertical gel unit will run 1-4 handcast gels in mini format, the sliding clamp assembly ensures fast set up times and leak-free operation. The insert for both gel casting and running eliminates time-consuming transfer of fragile gels. Designed for rapid turnaround of data, allowing PAGE, 2-D and blotting to be completed within a working day.

Technical Specifications	
Active gel size	8 x 8.5cm
Sample capacity	PAGE: 80 (20/Gel) Blot*: 4 2D*: 10 tubes
Buffer tank volume	Min: 250ml, Max: 1200ml
Unit dimensions (WxLxH)	19 x 13 x 15cm

Cat. No.	Code	Description
EP/70150	CVS10DSYS	OmniPAGE Mini, 10 x 10cm, 2 sets of glass plates with 1mm thick spacers, 2 x 12 sample combs, cooling pack, blanking plate and casting base
EP/70152	VS10DCAST	Dual caster, omniPAGE Mini, 10 x 10cm
EP/70154	CVS10EXCASTER	Casting stand, external without base
EP/70160	VS10NG	Glass plates, notched, 10 x 10cm, 2mm thick, pack of 2
EP/70162	VS10PG	Glass plates, plain, 10 x 10cm, 2mm thick, pack of 2
EP/70170	VS10-8MC-1	Comb, 8 sample, MC, 1mm thick
EP/70172	VS10-10-1	Comb, 10 sample, 1mm thick
EP/70174	VS10-12-1	Comb, 12 sample, 1mm thick
EP/70176	VS10-16MC-1	Comb, 16 sample, MC, 1mm thick
Accessories		
EP/70200	CVS1CBS	Blotting unit, complete system for mini vertical electrophoresis and blotting
EP/70202	VS10BI	OmniBlot mini insert including 4 x cassettes and 16 foam pads



EP/70150



EP/70160

*Will require purchase of additional items