

## Water bath RE400DB



A precise digital bath, specifically designed to be used with the Stuart rotary evaporators to give maximum efficiency for the smallest footprint. The clear LED digital display shows the actual bath temperature or the set temperature.

- Digital temperature display
- Corrosion resistant bowl
- Easy to clean
- Integrated carry handle
- Cool touch chemically resistant outer case

### Technical Specifications

Temperature range	Ambient to 90°C
Display resolution	1°C
Temperature control	±1°C
Capacity	4.3 Litres
Heater power	1400W
Dimensions	285 x 295 x 244mm

Cat. No.	Code	Description
ER/15030	RE400	Rotary evaporator complete with glass set 00
ER/15032	RE401	Rotary evaporator complete with glass set 01
ER/15034	RE402	Rotary evaporator complete with glass set 02
ER/15036	RE400P	Rotary evaporator complete with glass set 00P
ER/15038	RE401P	Rotary evaporator complete with glass set 01P
ER/15040	RE402P	Rotary evaporator complete with glass set 02P
ER/15024	RE400DB	Water bath, ambient to 90°C
ER/15026	SRC5	Recirculating cooler, -20°C
ER/15022	RE3022C	Rotary Vacuum pump, vacuum 9mbar
<b>Glassware and spares</b>		
ER/15106	FD50/4RE	Florentine flask 50ml (pack of 2)
ER/15108	FD100/4RE	Florentine flask 100ml
ER/15110	FD250/4RE	Florentine flask 250ml
ER/15112	FD500/4RE	Florentine flask 500ml
ER/15116	FD1L/4RE	Florentine flask 1litre
ER/15120	FD2L/4RE	Florentine flask 2 litre
ER/15124	RE100/RF/50	Receiving flask 50ml
ER/15126	RE100/RF/100	Receiving flask 100ml
ER/15128	RE100/RF/250	Receiving flask 250ml
ER/15130	RE100/RF/500	Receiving flask 500ml
ER/15132	RE100/RF/500P	Receiving flask 500ml plastic coated
ER/15134	RE100/RF/1L	Receiving flask 1litre
ER/15136	RE100/RF/1LP	Receiving flask, 1 litre, plastic coated
ER/15170	RE400/CO	Condenser, diagonal for glass set 00
ER/15172	RE400/COP	Condenser, diagonal plastic coated, for glass set 00P
ER/15174	RE400/VC	Condenser, vertical for glass set 01
ER/15176	RE400/VCP	Condenser, vertical, plastic coated, for glass set 01P
ER/15178	RE400/CF	Condenser, cold finger for glass set 02
ER/15180	RE400/CFP	Condenser, cold finger, plastic coated, for glass set 02P
ER/15182	RE400/VR	Feed tube/ vacuum release, all glass sets
ER/15154	RE400/VT	Vapour tube diagonal condensers, glass sets 00 & 00P
ER/15184	RE400/VT/CF	Vapour tube for vertical coil condensers, glass sets 01, 01P, 02, 02P
ER/15158	RE100VS	Vacuum seal, all glass sets
ER/15186	RE400/CONNECT	Connector, PTFE screw thread, all glass sets, pack of 10
ER/15188	RE400/ADPT	Adaptor for VC/CF condenser
JC35		Joint clip, spherical, for receiving flask, all glass sets (Pack of 2)
KC29		Joint clip, plastic for evaporating flask, all glass sets

Tubing for vacuum connection (8mm heavy wall) please order TU/76226

Tubing for condenser connection (8mm standard wall) please order TU/76108



ER/15034



ER/15034 & bath



ER/15024