

## Extraction glassware assemblies

PYREX®

Each assembly comprises the basic components for solid-liquid extraction.

- Flask
- Soxhlet extractor
- Associated lids and clips
- Condenser

Cat. No.	Code	Flask capacity, ml	Thimble size, mm
EX/12008	100RASX	250	30 x 100 (ET/17020)
EX/12010	200RASX	500	43 x 123 (ET/17024)

Large capacity Soxhlet extraction unit comprising of the following individual parts. Please note, for a full working assembly, all parts must be ordered.

Cat. No.	Description
EX5/75	Extractor Soxhlet 600ml, 75mm B34/35
DA5/75	Adapter multisocket lid 75mm
JC75F	Clip for 75mm flange
FR1L/5S	Flask round base short neck 1 litre 34/35
CX5/25/SC	Condenser Davies 200mm
ET/17026	Extraction thimble, single thickness 60 x 180



EX/12010

## Extraction thimbles, single thickness

For use in soxhlet extractors. 1 mm wall thickness.

Supplied in Pk Qty 25.



Cat. No.	Code	Size, dia.x length mm
ET/17002	2800-185	18 x 55
ET/17004	2800-199	19 x 90
ET/17006	2800-228	22 x 80
ET/17008	2800-258	25 x 80
ET/17010	2800-250	25 x 100
ET/17012	2800-266	26 x 60
ET/17014	2800-288	28 x 80
ET/17016	2800-280	28 x 100
ET/17018	2800-308	30 x 80
ET/17020	2800-300	30 x 100
ET/17022	2800-373	37 x 130
ET/17024	2800-432	43 x 123
ET/17026	2800-608	60 x 180
ET/17028	2800-226	22 x 65
ET/17030	2800-259	25 x 90
ET/17032	2800-260	26 x 100
ET/17034	2800-282	28 x 120
ET/17036	2800-330	33 x 100
ET/17038	2800-331	33 x 118
ET/17040	2800-338	33 x 80
ET/17042	2800-339	33 x 94
ET/17044	2800-412	41 x 123



ET/170