



FB/361

Filters, bottle top, PES

NALGENE®

Bottle top filters with 50, 75 or 90mm membranes. Securely screw onto media bottles with 33 or 45mm neck size. Allow direct filtration into sterile glass media bottles*. Moulded-in graduations. Side-arm features a quick-disconnect tubing adapter. Colour coded blue. Supplied sterile in cases of 12 units.

Cat. No.	Code	Capacity, ml	Membrane dia, mm	Pore size, µm
For 33mm neck bottles				
FB/36112	296-3345	150	50	0.45
FB/36114	295-3345	500	75	0.45
For 45mm neck bottles				
FB/36212	296-4545	150	50	0.45
FB/36214	295-4545	500	75	0.45

*Caution: Use only sterile bottle designed for use with vacuum. Always use safety shield during vacuum procedures.



FM/3300

Filter units, analytical, CN membrane

NALGENE®

Designed for ease of use this filter unit consists of a polypropylene upper and polystyrene lower with a white cellulose nitrate filter. Upper twists off easily for improved retrieval of the 47mm membrane. Wider base is stable on the bench. Non-silicone gasket reduces risk of contamination. Excellent for sterility testing, water quality or environmental testing. Supplied in sterile case of 12 units, cases of 72 are also available, please enquire.

Cat. No.	Code	Capacity, ml	Grid	Pore size, µm
FM/33002	130-4020	150	None	0.2
FM/33004	130-4045	150	Grey	0.45



FU/361

Filter funnel, analytical test, CN membrane

NALGENE®

Disposable, sterile filter funnels consisting of a polypropylene upper, high impact polystyrene collar and 47mm cellulose nitrate membrane. Designed for microbiological analysis testing. The grey membrane turns black when wet making it easier to detect and count growing colonies. Upper separates easily from lower for easy access to the membrane. Supplied sterile in cases of 50 with six funnel adapters in each case. Funnel adapter allows use with No. 8 stoppers.

Cat. No.	Code	Capacity, ml	Membrane	Pore size, µm	Pk Qty
FU/36132	145-0020	100	Plain	0.2	50
FU/36134	145-0045	100	Gridded	0.45	50
FU/36136	145-2020	250	Plain	0.2	50
FU/36138	145-2045	250	Gridded	0.45	50
FU/36140	147-0045	100	Grey, gridded	0.45	50
Accessories					
FU/36142	DS0396	Filter stopper			3
FU/36144	DS0397	Filter funnel adapter			25

Reusable Filterware

NALGENE®

Filter holder, autoclavable, polysulfone with graduated receiver, fitted with 2 side arms which accept 6 to 8mm i.d. vacuum tubing and also accepts syringe filters. The tubing adapter can be stuffed with cotton for sterile venting during autoclaving.



FU/362 reusable filters

Technical Specifications	300-4000	300-4050	300-4100
Upper chamber capacity, ml	250	500	500
Receiver capacity, ml	250	500	1000
Diameter of receiver, mm	83	117	114
Height including ports, mm	180	230	293
Width including side arms, mm	110	134	134

Cat. No.	Code	Description
FU/36232	300-4000	Filter holder and receiver 250 / 250ml
FU/36234	300-4050	Filter holder and receiver 500 / 500ml
FU/36240	300-4100	Filter holder and receiver 500 / 1000ml