



BA/40600



BA/40602



BA/40610

Balances, portable, Ohaus Scout SKX series



The Ohaus range of Scout SKX balances is built to endure the demanding workload of a busy laboratory environment. Combined with a simplified 4 button operation, the large backlit display increases readability in dim lighting conditions, and makes it easier to view results from across a desk or laboratory bench. Manufactured with ABS housing, stainless steel pan, ABS draft shield (1mg model only), weigh below hook, transportation lock, manual calibration lock switch and a security slot.

Connectivity options include RS232, USB Host, USB Device, Ethernet and Bluetooth®, (available as an extra), the Scout provides added versatility to your data collection and transfer needs.

- Stackable with storage cover (excluding 1 mg model)
- Stability indicator
- Mechanical and software overload/underload protection
- Low battery indicator
- Auto shut-off
- Auto tare
- Span calibration mass included on models up to 620g capacity
- AC adapter included (will run on 4 x AA batteries, not supplied)
- Expanded connectivity options including RS232, USB Host, USB Device, Ethernet and Bluetooth®
- Suitable for Density Determination and Mole Weighing

Cat. No.	Code	Description	Pan size
BA/40600	30253070	Balance, Ohaus Scout SKX123, 120g x 0.001g	93mm Dia.
BA/40602	30253071	Balance, Ohaus Scout SKX222, 220g x 0.01g	120mm Dia.
BA/40604	30253072	Balance, Ohaus Scout SKX422, 420 x 0.01g	120mm Dia.
BA/40606	30253073	Balance, Ohaus Scout SKX622, 620 x 0.01g	120mm Dia.
BA/40608	30253076	Balance, Ohaus Scout SKX421, 420 x 0.1g	120mm Dia.
BA/40610	30253077	Balance, Ohaus Scout SKX621, 620 x 0.1g	140 x 170mm
BA/40612	30253078	Balance, Ohaus Scout SKX2201, 2200 x .0.1g	140 x 170mm
BA/40614	30253079	Balance, Ohaus Scout SKX6201, 6200 x 0.1g	140 x 170mm
BA/40616	30253080	Balance, Ohaus Scout SKX8200, 8200 x 1g	140 x 170mm

Accessories

BA/40630	30268982	RS232 Interface Kit
BA/40632	30268983	USB Host Interface Kit
BA/40634	30268984	USB Device Interface Kit
BA/40636	30268985	Bluetooth Interface Kit
BA/40638	30268986	Ethernet Interface Kit
BA/40640	30268987	Stacking & Storage Cover (6 pcs)
BA/40642	30268988	Stacking & Storage Cover (1 pc)
BA/40646	30333960	Auxiliary display
BA/40650	30269020	Top Loading Kit for Density Determination
BA/40652	30269021	Carrying Case
BA/40654	30269022	In-use Cover
BA/40560	30064202	Impact printer SF40A
BA/40116	80850043	Security Device

Balances, portable, HCB Highland range



Adam Equipment's Highland balances have what it takes for lab work, field use and various education and industrial applications.

- HandiCal™ internal calibration for fast and accurate readings
- Internal rechargeable battery with auto power off and low-battery indicator
- ShockProtect™ overload protection
- RS-232 and USB interfaces
- Lock down mounting slot for Kensington™ type lock and cable
- Stackable for storage and transportation
- Removable draft shield 132mm diameter 100mm high
- Capacity tracker
- Audible overload indicator
- Pan diameter 120mm
- Dual tare keys, colour coded for easy recognition



BA/5042

Cat. No.	Model	Capacity	Readability
BA/50418	HCB 153	150g	0.005g
BA/50420	HCB 302	300g	0.01g
BA/50422	HCB 602H	600g	0.01g
BA/50424	HCB 602	600g	0.02g
BA/50426	HCB 1002	1000g	0.01g
BA/50428	HCB 1502	1500g	0.05g
BA/50430	HCB 3001	3000g	0.1g