



FC/62166

Flasks, conical, wide neck, borosilicate (Erlenmeyer)

Cat. No.	Capacity,ml	Neck I.D.,mm	Pk Qty
FC/62140	50	34	10
FC/62146	100	34	10
FC/62156	250	50	10
FC/62166	500	50	10
FC/62176	1000	50	10
FC/62180	2000	72	6



FC/124

Flasks, conical, wide neck, heavy duty

Reinforced rims for extra strength and resistance to impact.

PYREX®

Cat. No.	Code	Capacity,ml	Neck I.D.,mm	Pk Qty
FC/12448	5100-125	125	27	12
FC/12456	5100-250	250	33	12
FC/12466	5100-500	500	43	6
FC/12476	5100-1L	1000	48	6
FC/12480	5100-2L	2000	58	each



FC/126

Flasks, conical, screwcap

Ideal for storage of culture specimens and media. Phenolic screwcap with PTFE faced sealing disc. Flask and cap are autoclavable up to 121°C.

PYREX®

Cat. No.	Code	Capacity,ml	Cap size	Pk Qty
FC/12604	1137/13M	25	22	10
FC/12606	1137/02M	50	22	10
FC/12608	1137/03D	100	30	10
FC/12610	1137/06M	250	42	10
FC/12612	1137/08D	500	42	10
FC/12614	1137/10M	1000	42	each



FC/162

Flasks, conical, wide neck with screw closure

Polypropylene, or crystal clear PMP with PP screw cap. NS stoppers can also be used (not included). Flask and cap can be autoclaved at 121°C.

Azlon

Cat. No.	Code	Capacity,ml	Graduation intervals,ml	Neck size	Thread size
Polypropylene, translucent, blue graduations					
FC/16280	FDE542	50	10	34/35	GL40
FC/16282	FDE544	100	20	34/35	GL40
FC/16284	FDE546	250	50	45/40	GL52
FC/16286	FDE548	500	100	45/40	GL52
FC/16288	FDE550	1000	200	45/40	GL52
PMP, crystal clear, red graduations					
FC/16870	FDE442	50	10	34/35	GL40
FC/16872	FDE444	100	20	34/35	GL40
FC/16874	FDE446	250	50	45/40	GL52
FC/16876	FDE448	500	100	45/40	GL52
FC/16878	FDE450	1000	200	45/40	GL52