



FK/120

Flasks, Kjeldahl

Suitable for the determination of Nitrogen etc.

PYREX®

Cat. No.	Code	Capacity,ml	Pk Qty
FK/12046	1200/04M	100	10
FK/12066	1200/12D	500	10
FK/12074	1200/16M	800	10

Flasks, volumetric, Class B with stopper

Cat. No.	Capacity,ml
FV/52018	5
FV/52028	10
FV/52034	25
FV/52040	50
FV/52046	100
FV/52054	200
FV/52056	250
FV/52066	500
FV/52076	1000
FV/52080	2000
FV/52088	5000



FV/52066

Flasks, volumetric, Class B

Borosilicate glass with polyethylene stopper. Blue enamel inscriptions

MBL®

Cat. No.	Code	Capacity,ml	Stopper size	Tolerance ±ml	Pk Qty
FV/12018	FRD022P	5	10/19	0.05	5
FV/12028	FRD024P	10	10/19	0.05	5
FV/12032	FRD025P	20	10/19	0.08	5
FV/12034	FRD026P	25	10/19	0.08	5
FV/12040	FRD028P	50	12/21	0.12	5
FV/12046	FRD030P	100	14/23	0.2	5
FV/12054	FRD032P	200	14/23	0.3	5
FV/12056	FRD034P	250	14/23	0.3	2
FV/12066	FRD036P	500	14/23	0.5	2
FV/12076	FRD038P	1000	24/29	0.8	2
FV/12080	FRD040P	2000	29/32	1.2	each
FV/12088	FRD042P	5000	34/35	2.4	each



FV/120