

Balances, portable CQT Core series



The Core series from Adam Equipment is designed to withstand field and laboratory environments. Complete with an AC adapter, and removable draught shield on 120mm pan models.

- Three-point ShockProtect™ overload protection
- Built in security slot for Kensington™ type lock
- Large backlit display is clear and quick to read
- Stacking storage (CQT202, 251, 601 and 1501 models)
- Stainless steel pan
- Battery powered with 6 x AA batteries or AC adapter
- Removable draft shield on 120mm ø pans
- Auto power off to save battery life
- Low-battery indicator
- Dual tare keys, colour-coded for easy recognition



BA/5011

Cat. No.	Model	Capacity	Readability	Pan size
BA/50114	CQT202	200g	0.01g	120mm
BA/50116	CQT251	250g	0.1g	120mm
BA/50118	CQT601	600g	0.1g	120mm
BA/50120	CQT1501	1500g	0.1g	145mm
BA/50122	CQT2601	2600g	0.1g	145mm
BA/50124	CQT2000	2000g	1g	145mm
BA/50126	CQT5000	5000g	1g	145mm

Balances, Nimbus Precision Range



Designed for scientists, the Nimbus range of balances has a small footprint for economical use of laboratory bench space. The single piece construction offers greater stability producing highly repeatable results. The weighing sensor is crafted from a single block and contains fewer parts than a traditional force motor balance. This optimised internal mechanism is fabricated using hardened material, helping the Nimbus to tolerate rigorous laboratory use. Efficiency in the mechanics, plus improved processing power of the electronics, leads to excellent performance.

Suitable for weighing, parts counting, percentage weighing and density determination when used with the additional density kit.



BA/50854

Cat. No.	Code	Capacity	Readability	Pan diameter mm
BA/50854	NBL423e	420g	0.001g	120mm
BA/50856	NBL2602e	2600g	0.01g	160mm
BA/50858	NBL3602e	3600g	0.01g	160mm
BA/50860	NBL4602e	4600g	0.01g	160mm
BA/50862	NBL6201e	6200g	0.1g	160mm

Accessories

Cat. No.	Code	Description
BA/50866	1060012714	Density kit for 120mm and 160mm pan diameters for models up to 6Kg capacity
BA/50868	3012313008	Dust cover
BA/50872	3012313010	In use cover for 120mm pan
BA/50874	3012313011	In use cover for 160mm pan
BA/50876	1120011156	ATP Adam thermal printer
BA/50878	3126011263	Thermal printer roll
BA/50880	3014011014	RS-232 Cable
BA/50882	600002028	Adam DU data collection program
BA/50884	3074010267	USB Cable
BA/50886	3011413013	Weigh below hook for models up to 6Kg capacity



BA/50862