

## Funnels, filter, ribbed, long stem

**PYREX®**

With interior "ribs" to increase the filtration rate. Stem approximately 150mm long and 8mm OD

Cat. No.	Code	Dia.,mm	Pk Qty
FF/12306	2150/03M	65	10
FF/12308	2150/06M	80	10
FF/12310	2150/10M	110	10



FF/123

## Funnel, powder, borosilicate glass

**PYREX®**

With short wide stem for smooth flow.

Cat. No.	Code	Top dia., mm	Stem OD x length, mm	Pk Qty
FF/12408	2130/02M	80	18 x 25	10



FF/12408

## Funnels, general purpose, conical, HDPE

Tough, robust, general purpose funnels suitable for a wide range of liquids. Upstand on rim to reduce possible overflow. Ribbed on outside to eliminate back pressure.

Cat. No.	Top dia., mm	Pk Qty
FF/46102	80	5
FF/46104	100	5
FF/46106	120	5
FF/46108	140	5
FF/46110	180	5
FF/46112	260	5
FF/46114	420	5



FF/461

## Funnels, filter, PP

Body of funnel is exact 60° angle; with internal and external ribs.

Cat. No.	Code	Dia., mm	Stem, mm	Pk Qty
FF/46202	145	25	38	20
FF/46204	146	35	38	20
FF/46206	147	45	45	20
FF/46208	148	65	64	20
FF/46210	150	80	69	10
FF/46212	152	100	84	10
FF/46214	153	120	86	10
FF/46216	155	150	114	5
FF/46218	156	180	142	5



FF/4620