



GH/93002

## Gas cylinder supports, bench mounted

Designed to support gas cylinders up to 300 mm diameter. Heavy-duty alloy casting secured to benches up to 70 mm thick with two plated steel tommy screws. Two countersunk holes are provided to accept wood screws for permanent fixing, complete with fire resistant belt and buckle.

Cat. No.	Code	Description
GH/93002	3603	Bench mounted cylinder holder



GH/93004

## Gas cylinder supports, wall mounted

Comprised of a curved nylon coated steel bracket with a plated chain to hold cylinder firmly in position. For use with cylinders up to 230 mm diameter.

Cat. No.	Code	Description
GH/93004	9501	Wall mounted cylinder holder

## Gas Cylinder Trolleys

### Gas cylinder trolley, single

Designed to carry cylinders safely and to act as a stand to store the cylinder. The trolleys have 200mm steel centre wheels with rubber tyres. The larger trolley is designed to carry the larger LPG cylinders.

Cat. No.	Code	Height, mm	Width, mm	Depth, mm	Max Cylinder Diameter,mm
GT/14002	CYH100	1040	460	220	250
GT/14004	CYH101	1040	610	390	390



GT/14002

### Gas cylinder trolley, twin

Designed to move and store cylinders safely and with ease. The trolley has 355mm steel centre with rubber tyres.

Cat. No.	Code	Height, mm	Width, mm	Max Cylinder Diameter,mm
GT/14008	CYH201	1040	800	265/240



GT/14008