

## Haemocytometers


### Haemocytometer sets, to BS748

Complete set with counting chamber, cover glasses and one each red and white cell pipettes with tubing and mouthpiece in fitted plastic case. The metallized (bright-line) versions reduce glare and provide greater contrast between lines and background.

| Cat. No.   | Description             |
|--|-------------------------|
| <b>Improved Neubauer ruling, cell depth 0.1mm</b>      |                         |
| HA/13002   | Single cell, standard   |
| HA/13004   | Single cell, metallized |
| HA/13006   | Double cell, standard   |
| HA/13008   | Double cell, metallized |
| <b>Modied Fuchs-Rosenthal ruling, cell depth 0.2mm</b> |                         |
| HA/13010   | Double cell, metallized |
| <b>Thoma ruling, cell depth 0.1 mm</b>                 |                         |
| HA/13012   | Double cell, metallized |
| HA/13014   | Single cell, standard   |
| <b>Neubauer ruling</b>                                 |                         |
| HA/13016   | Single cell, standard   |
| HA/13018   | Single cell, metallized |
| HA/13020   | Double cell, standard   |
| HA/13022   | Double cell, metallized |
| <b>Burker ruling</b>                                   |                         |
| HA/13024   | Single cell, metallized |
| HA/13026   | Double cell, metallized |

#### Specialist counting chambers

| Cat. No. | Description   |
|----------|---|
| HA/13028 | Helber Bacteria 1 cell Thoma  |
| HA/13030 | Bacteria evaluation chamber   |
| HA/13032 | Malassez  |
| HA/13034 | Nageotte 2 cell   |
| HA/13036 | Spare cover glasses for Nageotte cell                                 |
| HA/13038 | Fertility Semen counting  |
| HA/13044 | Sedgewick rafter counting chamber, plastic, includes 1 x cover glass  |
| HA/13046 | Sedgewick rafter counting chamber, glass, includes 1 x cover glass    |
| HA/13052 | Sedgewick Rafter counting chamber without grid, including cover glass |
| HA/13048 | Cover glasses for Sedgewick rafter counting chamber, sold singly      |




HA/130

## Haemocytometers spares and accessories

### Cover glass

| Cat. No. | Description                                  | Pk Qty |
|----------|--|--------|
| HA/13122 | For single or double chamber sets. 22 x 25mm | 10     |



HA/1312

### Counting chambers

| Cat. No.   | Description                                | Pk Qty |
|--|--|--------|
| <b>Improved Neubauer ruling, 0.1mm cell depth</b>        |  |        |
| HA/13124   | Single cell, standard                      | 1      |
| HA/13126   | Single cell, metallized                    | 1      |
| HA/13128   | Double cell, standard                      | 1      |
| HA/13130   | Double cell, metallized                    | 1      |
| <b>Modified Fuchs-Rosenthal ruling, 0.2mm cell depth</b> |  |        |
| HA/13132   | Double cell, metallized                    | 1      |
| <b>Thoma ruling, 0.1mm cell depth</b>                    |  |        |
| HA/13134   | Double cell, metallized                    | 1      |
| <b>Parasitology chamber, Improved McMaster worm egg</b>  |  |        |
| HA/13136   | 10mm square with 10 equal vertical columns | 1      |
| HA/13146   | Hawkesley version                          |        |