



HY/65012

Hybridisation tubes/carriers and accessories for Hybridiser HB-1D & Hybrigene HB-3D

Adjustable speed rocking platform

Fits on the base of the unit still allowing hybridisation tubes to be placed above. Enables membranes to be processed using a wave motion. Available in fast or slow speeds, which range from 5 to 60 oscillations per minute.

Static shelf


Up to 2 (Hybrigene) or 3 (Hybridiser HB-1D) static stainless steel mesh shelves can be placed into the tube locators allowing the unit to be used as a simple non-motion incubator. Each shelf holds up to 4 kg.

Membrane separators

Strong, re-usable porous membranes for placing between blots. Ensure even hybridisation when processing multiple hybridisation membranes in one tube. Available as packs of 5 (20x20 cm).

Hybridisation glass tube holder

Useful storage facility for tubes during membrane loading or when not in use. Carries up to 3 large and 3 small tubes.



HY/65026

- Glass hybridisation tubes reduces volumes, simplifies washing, improves signals
- 3 sizes available: small, mini, and unique large tubes
- Durable: thick borosilicate glass
- Easy pour: non drip, safe and convenient
- Sealing Ring: "O" ring sealed end caps allow no leakage
- Broken thread on screw cap: easy opening with no vacuum problems
- Tube Assemblies: Adaptors to hold multiple small glass tubes to increase capacity

Cat. No.	Code	Description	Total capacity
Hybridiser HB-1D			
HY/65026	FHB12	Large glass tube	6
HY/65028	FHB11	Large tube, no cap	6
HY/65030	FHB32	Small glass tube	6
HY/65012	FHB41	Mini glass tube assembly	
HY/65034	FHB43	Tube carrier	6 x 4
HY/65036	F15ML1TH	15ml tube holder	6 x 16
HY/65038	F50ML1TH	50ml tube holder	4 x 8
HY/65056	FHB1/PLAT	Rocking platform 0-60 opm	1
HY/65050	FHB5H1	Static shelf, stainless steel mesh	1
Hybrigene HB-3D			
HY/65006	FHB16	Large glass tube with push-on end caps, 80 x 240mm	4
HY/65010	FHB36	Small glass tube with screw-on end caps, tube has no end plugs, 44 x 240mm	4
HY/65012	FHB41	Mini glass tube assembly with screw-on end cap, requires carrier FHB61	
HY/65014	FHB61	Tube carrier for mini glass tubes FHB41	4 x 4
HY/65022	F15ML4TH	Carrier with clips to hold 16 x 15ml tubes (approx diameter 16.5mm)	4 x 16
HY/65024	F50ML4TH	Carrier with clips to hold 8 x 50ml tubes (approx diameter 29.0mm)	4 x 8
HY/65058	FHB4/PLAT	Rocking platform 0-60 OPM	1
HY/65048	FHB4/SHI	Static shelf, Stainless Steel mesh	1
Common to both			
HY/65052	FMEM2020	Separation membranes pack of 5	
HY/65054	FHTRACK	Hybridisation glass tube holder	



HY/65036



HY/65054