

Jugs

Jugs, moulded graduations, high clarity PP

- Ergonomically designed jugs
- Integral thumb grip on handle allows easy carrying even when full and heavy.
- Made from tough durable polypropylene
- Autoclavable at 121°C.
- Clear moulded graduations
- Full size range for all requirements

Cat. No.	Code	Capacity, ml	Sub divisions, ml	Nominal dimensions (h x rim diam.), mm	Pk Qty
JU/16740	JPM0050P	50	10	60 x 50	5
JU/16746	JPM0100P	100	20	72 x 61	5
JU/16756	JPM0250P	250	50	100 x 90	1
JU/16766	JPM0500P	500	50	120 x 102	1
JU/16776	JPM1000P	1000	100	180 x 122	1
JU/16780	JPM2000P	2000	100	189 x 151	1
JU/16784	JPM3000P	3000	200	215 x 172	1
JU/16788	JPM5000P	5000	500	253 x 199	1



JU/167

Jugs, moulded graduations, short form, PP

Autoclavable, highly translucent especially when filled with liquids. Permanently moulded in graduations. Excellent chemical resistance. Non drip spout. Thumb grip on handle.

Cat. No.	Code	Capacity, ml	Graduations, ml
JU/46266	1156	500	25
JU/46276	1157	1000	50
JU/46280	1158	2000	50
JU/46284	1159	3000	100
JU/46288	1160	5000	250



JU/462

Jugs, printed graduations, high clarity PP

Azlon

- Printed jugs in high clarity polypropylene
- High clarity PP is autoclavable at 121°C
- Ergonomic handle means comfortable safe handling of all sizes.
- Sharp blue graduations are easy to read and are indelible.
- Precise non-drip pouring spouts reduce splashing.
- Complete size range for all requirements

Cat. No.	Code	Capacity, ml	Sub divisions, ml	Nominal dimensions (h x rim diam.), mm	Pk Qty
JU/16240	JPB0050P	50	10	60 x 50	5
JU/16246	JPB0100P	100	20	72 x 61	5
JU/16256	JPB0250P	250	50	90 x 100	1
JU/16266	JPB0500P	500	50	120 x 102	1
JU/16276	JPB1000P	1000	100	150 x 122	1
JU/16280	JPB2000P	2000	100	187 x 151	1
JU/16284	JPB3000P	3000	200	215 x 172	1
JU/16288	JPB5000P	5000	500	253 x 199	1



JU/162

Jugs, stainless steel, graduated, tapered with spout

Cat. No.	Capacity, ml
JU/94066	500
JU/94076	1000
JU/94080	2000



JU/940