



Tape, labelling, cryogenic to -80°C

Laboratory tape, for labelling and easy identification of samples. Available in 4 sizes and 6 different colours. The paper based surface is suitable for all types of marker - fibre tips, ball point and gel pens. Ideal for use on most smooth plastic, glass and metal surfaces and can withstand temperatures down to -80°C.

Cat. No.	Description
LA/10302	Tape, 12m x 12.7mm, cryogenic, White, core size 2.7cm
LA/10304	Tape, 12m x 12.7mm, cryogenic, Red, core size 2.7cm
LA/10306	Tape, 12m x 12.7mm, cryogenic, Blue, core size 2.7cm
LA/10308	Tape, 12m x 12.7mm, cryogenic, Green, core size 2.7cm
LA/10310	Tape, 12m x 12.7mm, cryogenic, Yellow, core size 2.7cm
LA/10312	Tape, 12m x 12.7mm, cryogenic, Orange, core size 2.7cm
LA/10314	Tape, 12m x 19mm, cryogenic, White, core size 2.7cm
LA/10316	Tape, 12m x 19mm, cryogenic, Red, core size 2.7cm
LA/10318	Tape, 12m x 19mm, cryogenic, Blue, core size 2.7cm
LA/10320	Tape, 12m x 19mm, cryogenic, Green, core size 2.7cm
LA/10322	Tape, 12m x 19mm, cryogenic, Yellow, core size 2.7cm
LA/10324	Tape, 12m x 19mm, cryogenic, Orange, core size 2.7cm
LA/10326	Tape, 55m x 12.7mm, cryogenic, White, core size 7.7cm
LA/10328	Tape, 55m x 12.7mm, cryogenic, Red, core size 7.7cm
LA/10330	Tape, 55m x 12.7mm, cryogenic, Blue, core size 7.7cm
LA/10332	Tape, 55m x 12.7mm, cryogenic, Green, core size 7.7cm
LA/10334	Tape, 55m x 12.7mm, cryogenic, Yellow, core size 7.7cm
LA/10336	Tape, 55m x 12.7mm, cryogenic, Orange, core size 7.7cm
LA/10338	Tape, 55m x 19mm, cryogenic, White, core size 7.7cm
LA/10340	Tape, 55m x 19mm, cryogenic, Red, core size 7.7cm
LA/10342	Tape, 55m x 19mm, cryogenic, Blue, core size 7.7cm
LA/10344	Tape, 55m x 19mm, cryogenic, Green, core size 7.7cm
LA/10346	Tape, 55m x 19mm, cryogenic, Yellow, core size 7.7cm
LA/10348	Tape, 55m x 19mm, cryogenic, Orange, core size 7.7cm



LA/10316



LA/10346

Dispensers for labelling tape reels

Tape dispensers for labelling tape, available in 2 sizes, small for reels with 2.7cm core and a large version for tapes with a 7.7cm core. Both units are of solid metal construction and feature suction cups on the bottom of the unit to securely hold the dispenser to the benchtop. Both dispensers feature 4 positions that are easy to load and remove the tape reels and provide access to high volumes of tape from one handy location.

Cat. No.	Description
LA/10360	Dispenser, tape, small, 4 place, for tapes with 2.7mm core
LA/10362	Dispenser, tape, large, 4 place, for tapes with 7.7mm core



LA/10360 and LA/10362

Notes
