

Replacement filters:


All models require 1 main filter. Models VLFT1800, HLT1500 and HLFT1800 require 2 prefilters, all other models require one prefilter only.

Filters:

Cat. No.	Code	Description
LA/70060	K-PF0022	Pre-filter, replacement for laminar flow cabinets
LA/70062	K-HF0135	Filter, main UPLA, replacement for model VLFT1000
LA/70064	K-HF0136	Filter, main UPLA, replacement for model VLFT1200
LA/70066	K-HF0137	Filter, main UPLA, replacement for model VLFT1500
LA/70068	K-HF0138	Filter, main UPLA, replacement for model VLFT1800
LA/70070	K-HF0140	Filter, main UPLA, replacement for model HLFT1000
LA/70072	K-HF0141	Filter, main UPLA, replacement for model HLFT1200
LA/70074	K-HF0142	Filter, main UPLA, replacement for model HLFT1500
LA/70076	K-HF0143	Filter, main UPLA, replacement for model HLFT1800

Additional options:

Cat. No.	Code	Description
LA/70080	K-VLF1000ST	Base stand for model VLFT1000 (A Stainless Steel work surface must also be purchased)
LA/70082	K-VLF1200ST	Base stand for model VLFT1200 (A Stainless Steel work surface must also be purchased)
LA/70084	K-VLF1500ST	Base stand for model VLFT1500 (A Stainless Steel work surface must also be purchased)
LA/70086	K-VLF1800ST	Base stand for model VLFT1800 (A Stainless Steel work surface must also be purchased)
LA/70090	K-HLF1000ST	Base stand for model HLFT1000
LA/70092	K-HLF1200ST	Base stand for model HLFT1200
LA/70094	K-HLF1500ST	Base stand for model HLFT1500
LA/70096	K-HLF1800ST	Base stand for model HLFT1800
LA/70100	K-VLF1000WS	Work surface, Stainless steel for model VLFT1000
LA/70102	K-VLF1200WS	Work surface, Stainless steel for model VLFT1200
LA/70104	K-VLF1500WS	Work surface, Stainless steel for model VLFT1500
LA/70106	K-VLF1800WS	Work surface, Stainless steel for model VLFT1800
LA/70110	K-VLF1000UV	UV light, night door and stainless steel work surface for model VLFT1000
LA/70112	K-VLF1200UV	UV light, night door and stainless steel work surface for model VLFT1200
LA/70114	K-VLF1500UV	UV light, night door and stainless steel work surface for model VLFT1500
LA/70116	K-VLF1800UV	UV light, night door and stainless steel work surface for model VLFT1800
LA/70120	K-CDES	Electrical socket, double for models HLFT1000, HLFT1200, HLFT1500, HLFT1800
LA/70122	K-HLF1000UV	UV light and night door for model HLFT1000
LA/70124	K-HLF1200UV	UV light and night door for model HLFT1200
LA/70126	K-HLF1500UV	UV light and night door for model HLFT1500
LA/70128	K-HLF1800UV	UV light and night door for model HLFT1800



LA/70026