

Anti-vibration table

Providing a solid, stable surface developed specifically to help decrease vibration and provide a marked precision despite air currents or movements during weighing processes. Supplied for easy self-assembly, with all necessary tools provided, the unit comes in 2 parts, the table style base, and a scratch-resistant dark granite slab with a working surface measuring 400 x 450mm. The work surface sits on shock absorbing rubber mounts fitted to the base which also has levelling feet to correct uneven surfaces.

Cat. No.	Description
BA/55002	Anti-vibration table with granite work surface



BA/55002

Calibration Weights

Calibration weights, available individually or in sets to either Class M1 or F1 tolerance in a plastic case. Individual weights are supplied with COPRAC (Comite Francais d'Accreditation) certification. Sets are available with or without certification. Weights from 1mg to 500mg are manufactured from aluminium with weights from 1g to 500g being manufactured from polished stainless steel.

Please note, prices for E1, E2 and F2 tolerances are available on request

Individual weights:

Individual weights: Mass	M1 Cat. No.	F1 Cat. No.
1mg	BW/84100	BW/84200
2mg	BW/84102	BW/84202
5mg	BW/84104	BW/84204
10mg	BW/84106	BW/84206
20mg	BW/84108	BW/84208
50mg	BW/84110	BW/84210
100mg	BW/84112	BW/84212
200mg	BW/84114	BW/84214
500mg	BW/84116	BW/84216
1g	BW/84118	BW/84218
2g	BW/84120	BW/84220
5g	BW/84122	BW/84222
10g	BW/84124	BW/84224
20g	BW/84126	BW/84226
50g	BW/84128	BW/84228
100g	BW/84130	BW/84230
200g	BW/84132	BW/84222
500g	BW/84134	BW/84234

Uncertified weight sets:

Cat. No.	Description
BW/84300	Calibration weight set, M1 set of 12 range 1mg to 500mg
BW/84302	Calibration weight set, M1 set of 23 range 1mg to 200g
BW/84304	Calibration weight set, M1 set of 11 range 1g to 200g
BW/84306	Calibration weight set, M1 set of 13 range 1g to 1Kg
BW/84308	Calibration weight set, F1 set of 12 range 1mg to 500mg
BW/84310	Calibration weight set, F1 set of 23 range 1mg to 200g
BW/84312	Calibration weight set, F1 set of 11 range 1g to 200g
BW/84314	Calibration weight set, F1 set of 13 range 1g to 1Kg

COPRAC Certified weight sets:

Cat. No.	Description
BW/84320	Calibration weight set, M1 certified, set of 12 range 1mg to 500mg
BW/84322	Calibration weight set, M1 certified set of 23 range 1mg to 200g
BW/84324	Calibration weight set, M1 certified set of 11 range 1g to 200g
BW/84326	Calibration weight set, M1 certified set of 13 range 1g to 1Kg
BW/84328	Calibration weight set, F1 certified set of 12 range 1mg to 500mg
BW/84330	Calibration weight set, F1 certified set of 23 range 1mg to 200g
BW/84332	Calibration weight set, F1 certified set of 11 range 1g to 200g
BW/84334	Calibration weight set, F1 certified set of 13 range 1g to 1Kg



BW/84100



BW/84300



BW/84320