




MS/99002

Microscope slide storage and handling

Microscope slide box, wooden

Wooden box with hinged lid.

Cat. No.	Description
MS/99002	Box for 50 slides
MS/99004	Box for 100 slides



MS/9600

Microscope slide box, plastic

Designed for easy handling of 75 x 25mm slides. The 25 place slide boxes are two-piece and stack securely on top of each other. The 100 place slide boxes have clearly numbered slots which correspond to the index card in the lid. Secured by a rust resistant clasp these hinged boxes also stack securely on top of each other.

Cat. No.	No. places	Dimensions (wxdxh), mm
MS/96002	25	90 x 140 x 35
MS/96004	100	210 x 167 x 35

Microscope slide folder

Made from strong card with hinged flaps and numbered positions. Compartments are recessed to protect slides when the cover is closed. Holds 20 slides 76 x 26mm

Cat. No.	Description
MS/09680	Microscope slide folder




MS/96006

Microscope slide trays, cardboard

Cardboard microscope slide trays to hold 30 standard 75x25mm. The horizontal recessed compartments, which are thumb cut for easy removal of slides, protect specimens when the cover is closed. Pack of 40 trays.

Cat. No.	Description
MS/96006	Cardboard microscope slide tray to hold 30 standard slides, Black



MS/43002

Microscope slide dispenser

For finger print free dispensing of 76x26mm microscope slides. Up to 60 slides can be put in the magazine. Each slide is dispensed by rotating the knurled knob on either side of the moulded base.

Cat. No.	Code	Dimensions (wxdxh), mm
MS/43002	540	100 x 120 x 140

Kartell