



MS/16202

Microscope slide staining jar

Robust polypropylene jars, can hold up to 10 standard 76 x 26mm slides back to back.

Can also be used as a mailing container. Supplied complete with lid.

Azlon

Cat. No.	Code	Description
MS/16202	SWM012	Staining jar, square pattern



MS/92116

Microscope slide staining trough

Glass trough with loose flat lid 85 x 105 x 75mm high. Accepts one glass accessory tray which has 10 inner grooves to accept 2 slides 76 x 26mm horizontally back to back. Must be ordered separately, capacity 20 slides.

Cat. No.	Description
MS/92116	Staining trough
MS/92118	Accessory tray



MS/92204

Microscope slide staining trough, Hellendahl

Glass trough with loose flat lid. Internally grooved to accept 16 slides 76 x 26mm vertically back to back.

Cat. No.	Description
MS/92204	Staining trough

Microscope slide staining trough, Schiefferdecker

Glass trough with loose flanged lid. Internally grooved to accept 10 slides 76 x 26mm horizontally back to back, capacity 20 slides.

Cat. No.	Description
MS/92206	Staining trough



MS/26002

Blotting/drying sheets

Blocks of 50 sheets of absorbent paper 37 x 100mm strips for the removal of excess liquids, e.g. from microscopic specimens.



Cat. No.	Size, mm	Pk Qty
MS/26002	37 x 100	50