



MI/75202 with MI/75206



MI/7521

## Microfine grinder drive, MF 10 basic

Continuously operating universal grinder.

- Powerful drive
- Easy to clean working surface made of stainless steel
- Two different grinding heads can be attached to the drive
- Heads are easily changeable
- Grinding heads not incl. with delivery

Cat. No.	Code	Description
MI/75202	MF10	Microfine grinder drive
<b>Accessories</b>		
MI/75204	MF 10.1	Cutting grinding head
MI/75206	MF 10.2	Impact grinding head

## Sieves for MF 10 basic

Interchangeable sieves for insertion into the grinding heads ensure max. particle size filtering.

Cat. No.	Code	Description
MI/75212	MF 0.25	Sieve with 0.25mm holes
MI/75214	MF 0.5	Sieve with 0.50mm holes
MI/75216	MF 1.0	Sieve with 1.00mm holes
MI/75218	MF 2.0	Sieve with 2.00mm holes
MI/75220	MF 3.0	Sieve with 3.00mm holes

## Notes