



MB/85108

Homogeniser bags

BagLight® Plain Blender Bag, Sterile

BagLight® is a plain bag used for products not requiring filtration. Made from PolySilk®, it is remarkably resistant and can be deep-frozen.

Cat. No.	Code	Description	Dimensions	Pk Qty
MB/85108	132 025	BagLight® 400	180 x 290mm	500
MB/85110	133 025	BagLight® 3500	380 x 510mm	250

BagFilter® Side Filter Bag, Sterile

BagFilter® has a non-woven filter and is recommended for the analysis of raw and cooked meat, fish, vegetables, fruit etc. It has excellent mechanical resistance. After blending, the fibres are retained in one of the compartments leaving the filtered solution easily accessible for pipetting. Each case contains 20 sterile packs of 25.

Cat. No.	Code	Description	Dimensions	Pk Qty
MB/85112	111 425	BagFilter® 400	190 x 300mm	500
MB/85114	252 019	Bagtips straws		1000



MX/15326



MX/15328



MX/15330

Mixers / rollers

Mixers, roller, rocking action



These roller mixers provide a gentle, but highly efficient, rocking and rolling action. There are four models available, two analogue models with either 6 or 9 rollers, with a fixed speed of 33 rpm. Or two more advanced digital models with either 6 or 9 rollers, with digitally controlled speed between 5 and 60 rpm and programmable timer.

Recommended for mixing blood samples, viscous substances and liquid-solid suspensions where minimum aeration is required. Units can be used in incubators up to 60°C and humidity up to 80%, or in cold rooms down to 4°C.

- Rocking and rolling action for complete mixing
- Choice of analogue fixed speed or digital variable speed model
- Six or nine roller design with footprint
- Designed for continuous quiet operation
- With BioCote® antimicrobial protection

Technical Specifications

	Analogue	Digital
Number of rollers	6 (SRT6) 9 (SRT9)	6 (SRT6D) 9 (SRT9D)
Speed	33rpm	5 to 60rpm
Amplitude	16mm	16mm
Maximum load	5kg	5kg
Controls	Analogue	Digital
Timer	No	Yes
Roller size (l x d)	340x30mm	340x30mm
Overall dimensions	565x240x110mm (SRT6) 565x360x110mm (SRT9)	565x240x110mm (SRT6D) 565x360x110mm (SRT9D)
Net weight	5.1kg (SRT6) 6.9kg (SRT9)	5.1kg (SRT6D) 6.9kg (SRT9D)
Electrical Supply	230V, 50Hz, 50W	230V, 50Hz, 50W

Cat. No.	Code	Description
MX/15324	SRT6	Roller mixer, 6 rollers, analogue control, fixed speed
MX/15326	SRT6D	Roller mixer, 6 rollers, digital control, variable speed
MX/15328	SRT9	Roller mixer, 9 rollers, analogue control, fixed speed
MX/15330	SRT9D	Roller mixer, 9 rollers, digital control, variable speed