

Mortars and pestles

Mortars and pestles, glass

Dimension is the nominal outside diameter of the mortar.

| Cat. No. | Diameter |
|----------|----------|
| MP/92108 | 75mm |
| MP/92110 | 100mm |
| MP/92112 | 150mm |



MP/92108

Mortar and pestle, Agate

| Cat. No. | Diameter O.D | Capacity |
|----------|--------------|----------|
| MP/98144 | 50mm | 10ml |
| MP/98148 | 70mm | 25ml |
| MP/98152 | 95mm | 45ml |



MP/98144

Mortar and pestle, porcelain

With unglazed grinding surfaces and glazed external surfaces. Dimension is the nominal outside diameter of the mortar.

| Cat. No. | Capacity, ml | Diameter, mm |
|----------|--------------|--------------|
| MP/98126 | 30 | 70 |
| MP/98128 | 75 | 91 |
| MP/98130 | 150 | 113 |
| MP/98132 | 220 | 125 |
| MP/98134 | 325 | 144 |
| MP/98136 | 1135 | 215 |



MP/98130

Mortar and pestle, porcelain

Unglazed. Dimension is the nominal outside diameter of the mortar. Spare pestles available.

| Cat. No. | Diameter, mm |
|----------|--------------|
| MP/98120 | 115 |
| MP/98122 | 135 |
| MP/98124 | 150 |



MP/9812

Pestle Mixer, motorised

Pestle Mixer, lightweight battery powered for use with disposable pestles, mixers and micro tubes. Cell tissue and pellets are completely homogenized within seconds by means of a high speed vortexing action. Greatly reduces hand fatigue, comes complete with 2 "AA" batteries and Pestle/Mixer adapter.

Pestles & Tubes

Designed to homogenize cells and plant tissue or to re-suspend protein and DNA pellets. Individually wrapped to guarantee purity, DNase/RNase and pyrogen free. The grooved grip handle is easy to turn and does not slip even when wearing gloves. Autoclavable, pestle length 7.5cm.

| Cat. No. | Description |
|----------|--|
| MP/10000 | Pestle mixer, motorised, battery operated, for micro tubes |
| MP/10002 | Adapter, rubber, replacement, for connecting pestle to mixer |
| MP/10004 | Pestle, micro, 0.5ml |
| MP/10006 | Pestle, micro, 1.5ml |
| MP/10008 | Microtube, 1.5ml |
| MP/10010 | Microtube, 0.5ml |
| MP/10012 | Microtube and pestle together in one sleeve, 1.5ml |
| MP/10014 | Microtube and pestle together in one sleeve, 0.5ml |
| BA/93006 | Battery, AA, spare for Pestle Mixer, 2 required |



MP/10000