

Weighing bottles



Glass, tall form with hollow blown stopper. Pack of 4.

Cat. No.	Code	Capacity, ml	HT x Dia. mm
BB/32202	L180/B/4	15	50 x 25
BB/32204	L180/C/4	28	60 x 30
BB/32206	L180/D/4	60	80 x 40
BB/32208	L180/E/4	135	100 x 50



BB/3220

Weighing scoops, borosilicate glass



- Flat base for stability on the balance pan
- Hollow stem for pouring powders into flasks etc
- Ideal for weighing small quantities of powder

Cat. No.	Code	Capacity, ml	Pk Qty.
BS/12202	2115/02M	3	10
BS/12204	2115/04M	6	10
BS/12206	2115/06M	10	10



BS/1220

Basins

Basins, borosilicate glass, crystallising, flat bottom, without spout



Cat. No.	Code	Capacity, ml	Dia x H., mm	Pk Qty.
BC/12002	1470/02D	100	70x40	10
BC/12004	1470/04D	150	80x45	10
BC/12006	1470/08D	300	95x55	10
BC/12008	1470/10D	500	115x65	10
BC/12010	1470/12D	900	140x75	10
BC/12012	1470/16D	2000	190x90	10



BC/120

Basins, borosilicate glass, crystallising, flat bottom, without spout

Cat. No.	Capacity, ml	Dia, mm	Pk Qty.
BC/62102	40	50	10
BC/62104	60	60	10
BC/62106	100	70	10
BC/62108	150	80	10
BC/62110	300	95	10
BC/62112	500	115	10
BC/62114	900	140	10
BC/62116	2000	190	4



BC/62110

Basins, borosilicate glass, crystallising, flat bottom, with spout

Cat. No.	Capacity, ml	Dia, mm	Pk Qty.
BC/62002	40	50	10
BC/62004	60	60	10
BC/62006	100	70	10
BC/62008	150	80	10
BC/62010	300	95	10
BC/62012	500	115	10
BC/62014	900	140	10
BC/62016	2000	190	4



BC/62012