

Ovens, Apex and Peak Series (continued)

Technical Specifications	Peak Ovens			
Natural circulation, models				
	PN30	PN60	PN120	PN200
Volume (litres)	27	57	115	215
Temperature range	Ambient + 30 to 300°C			
Temperature stability (°C)	±0.5			
Temperature uniformity (°C)	±7.0@300°C			
Heat up time (mins)	52	52	52	58
Recovery time (mins)	8.5	8.5	8.5	10
Shelves fitted	2	2	2	2
Number of shelf positions	3	5	9	15
Dimensions:				
Internal HxWxD (mm)	255x330x320	350x392x420	450x492x520	700x592x520
External HxWxD (mm)	470x665x470	570x765x570	670x865x670	920x965x670
Weight (Kg)	30	45	60	75
Power supply	220V - 240V, 50 - 60Hz, single phase			
Max Power (W)	750	1000	1500	2250



OV/35004



OV/35012

Technical Specifications	Peak Ovens			
Fan circulated, model number:				
	PF30	PF60	PF120	PF200
Volume (litres)	28	66	127	230
Temperature range	Ambient + 30 to 300°C			
Temperature stability (°C)	±0.2			
Temperature uniformity (°C)	±5.0@300°C			
Heat up time (mins)	40	36	35	42
Recovery time (mins)	4	4	4	5
Air changes	10	21	11	10
Air changes when fitted with optional exhaust fan	20	153	79	50
Shelves fitted	2	2	2	2
Number of shelf positions	3	5	9	15
Dimensions:				
Internal HxWxD (mm)	300x290x320	400x390x420	500x490x520	750x590x520
External HxWxD (mm)	470x665x470	570x765x570	670x865x670	920x965x670
Weight (Kg)	30	45	60	75
Power supply	220V - 240V, 50 - 60Hz, single phase			
Max Power (W)	1000	1500	2000	2700



OV/35016

Cat. No.	Code	Description
OV/35002	PN030	Oven, 27 L
OV/35004	PF030	Oven with fan, 28 L
OV/35006	PN060	Oven, 57 L
OV/35008	PF060	Oven with fan, 66 L
OV/35010	PN0120	Oven, 115 L
OV/35012	PF0120	Oven with fan, 127 L
OV/35014	PN0200	Oven, 215 L
OV/35016	PF0200	Oven with fan, 230 L

Notes
