



BE/120

Basins, borosilicate glass, evaporating, flat bottom, with spout

**PYREX®**

Cat. No.	Code	Capacity, ml	Dia x H, mm	Pk Qty
BE/12002	1450/01D	15	50x25	10
BE/12004	1450/03D	45	60x30	10
BE/12006	1450/04D	90	80x45	10
BE/12008	1450/06D	170	95x55	10
BE/12010	1450/10D	320	115x65	10
BE/12012	1450/14D	600	140x80	10



BE/622

Basins, borosilicate glass, evaporating, flat bottom with spout

Cat. No.	Capacity, ml	Dia, mm	Pk Qty
BE/62202	45	60	10
BE/62204	90	80	10
BE/62206	170	95	10
BE/62208	320	115	10



BE/932

Basin, glazed porcelain, evaporating, round bottom, shallow, with spout

Cat. No.	Capacity, ml	Dia, mm	Height, mm	Pk Qty
BE/93202	41	72	26	5
BE/93204	62	83	27	5
BE/93206	110	96	30	5
BE/93210	215	117	41	1
BE/93212	240	127	42	1
BE/93214	325	142	46	1
BE/93216	420	154	50	1
BE/93218	558	170	55	1



BE/933

Basin, glazed porcelain, evaporating, round bottom, deep, with spout

Cat. No.	Capacity, ml	Dia, mm	Height, mm	Pk Qty
BE/93302	36	61	28	5
BE/93304	72	77	33	5
BE/93306	107	81	38	5
BE/93308	140	94	42	1
BE/93310	190	103	45	1
BE/93312	232	112	47	1
BE/93316	500	138	61	1
BE/93318	590	148	66	1
BE/93320	920	170	74	1
BE/93322	1135	180	84	1
BE/93324	1500	202	86	1



BE/935

Basin, glazed porcelain, evaporating, flat bottom, with spout

Cat. No.	Capacity, ml	Dia, mm	Height, mm	Pk Qty
BE/93502	22	54	22	5
BE/93504	62	70	30	5
BE/93506	93	86	33	5
BE/93508	154	98	40	10
BE/93510	265	112	50	10
BE/93512	330	126	53	10
BE/93514	450	145	55	10
BE/93516	620	151	63	10
BE/93518	800	169	68	10
BE/93520	1200	195	74	1
BE/93522	3000	254	105	1