

## Pipettes, graduated

Pipettes, Graduated, Class AS, soda glass,  
to BS700, ISO 835



Amber stain graduations, individually calibrated and printed. Type 3 - Calibrated to deliver from zero at the top down to the jet, run-out.

Cat. No.	Capacity, ml	Graduation ml	Pk Qty
PG/32310	1	0.01	5
PG/32312	2	0.02	5
PG/32318	5	0.05	5
PG/32328	10	0.1	5
PG/32334	25	0.2	5



PG/32

## Pipettes, graduated, class B

**MBL®**

- Complies with ISO 835
- Graduations in permanent amber stain which withstands aggressive cleaning solutions
- Manufactured from soda-lime glass with robust tooled jets

Type 1 - Calibrated to deliver from zero at the top to any graduation line down to the shoulder

Type 2 - Calibrated to deliver from any graduation line down to zero at the jet

Type 3 - Calibrated to deliver from zero at the top to any graduation line down to the jet

Type 4 - Calibrated to deliver from any graduation line down to zero at the jet with the last drop expelled by blowing

Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cap., ml	Grad., ml	Pk Qty
Type 1	Type 2	Type 3	Type 4			
PG/12110	PG/12210	PG/12310	PG/12410	1	0.01	5
PG/12112	PG/12212	PG/12312	PG/12412	2	0.02	5
PG/12118	PG/12218	PG/12318	PG/12418	5	0.05	5
PG/12128	PG/12228	PG/12328	PG/12428	10	0.1	5
PG/12134	PG/12234	PG/12334	PG/12434	25	0.2	5



PG/12128



PG/122

## Pipettes, graduated, class AS

**MBL®**

- Complies with ISO 835
- Graduations in permanent amber stain which resists aggressive washing solutions
- Calibrated 'to deliver' with a +5s waiting time
- Manufactured from soda-lime glass with tooled jets for extra strength and to ensure correct outflow times

Type 1 - Calibrated to deliver from zero at the top to any graduation line down to the shoulder

Type 2 - Calibrated to deliver from any graduation line down to zero at the jet

Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cap., ml	Grad., ml	Tol., ± ml	Pk Qty
Type 1	Type 2	Type 1 with cert*	Type 2 with cert*				
PG/12510	PG/12610	PG/12710	PG/12810	1	0.01	0.007	2
PG/12512	PG/12612	PG/12712	PG/12812	2	0.02	0.01	2
PG/12518	PG/12618	PG/12718	PG/12818	5	0.05	0.03	2
PG/12528	PG/12628	PG/12728	PG/12828	10	0.1	0.05	2
PG/12534	PG/12634	PG/12734	PG/12834	25	0.2	0.1	2

\*Supplied with a traceable calibration certificate at 5 points



PG/12528



PG/128

## Pipettes, Dropping, Neutral glass, with Teat

Cat. No.	Description	Pk Qty
PP/55000	Dropping pipette with teat graduated at 0.5 and 1.0ml	10
PP/55002	Dropping pipette with teat 2ml	10
PP/55004	Dropping pipette without teat graduated at 0.5 and 1.0ml	10
PP/55006	Teat, black natural rubber	10



PP/5500