



PS/16



PS/165



PS/96



PP/12040

Pipettes, Serological

Pipettes, serological, disposable, polystyrene



Sterilin pipettes are available either bulk or individually packaged. The pipettes feature both negative and reverse graduations and have a non-wettable polystyrene surface to prevent readings from being distorted by a meniscus effect. Certified non-pyrogenic to <math><0.25\text{EU/ml}</math> with an accuracy of $\pm 1\%$ the pipettes are assured sterile by gamma irradiation.

Supplied with suction adaptor unless stated otherwise.

Pipettes Standard – Plastic Peel and Bulk

Cat. No.	Code	Description	Inner Pk	Pk Qty
PS/16110	40101K	1 ml pipette, single wrap (shelf pack 100)	1	1000
PS/16112	40102K	2 ml pipette, single wrap (shelf pack 100)	1	1000
PS/16118	40105	5 ml pipette, single wrap (shelf pack 100), no adapter	1	200
PS/16218	47105	5 ml pipette, single wrap with suction adaptor (shelf pack 50)	1	200
PS/16134	40125	25 ml pipette, single wrap with suction adaptor (shelf pack 25)	1	200
PS/16140	47150	50 ml pipette, single wrap with suction adaptor (shelf pack 50)	1	50
PS/16310	40301K	1 ml pipette	50	1000
PS/16312	40302K	2 ml pipette	50	1000
PS/16318	40305	5 ml pipette, no adapter	25	200
PS/16428	47310	10 ml pipette, no adapter	25	500

Pipettes Standard - Paper Peel



Paper peel packaging printed with product code, lot number and expiry date.

Cat. No.	Code	Description	Inner Pk	Pk Qty
PS/16510	40501K	1ml pipette, single wrap	50	1000
PS/16512	40502K	2ml pipette, single wrap	50	1000
PS/16518	47505	5ml pipette, single wrap	50	200
PS/16528	47510	10ml pipette, single wrap	50	500
PS/16534	47525	25ml pipette, single wrap	25	200
PS/16540	47550	50ml pipette, single wrap	25	50

Pipettes Serological, Disposable, Polystyrene

Single use made of crystal polystyrene, sterile.

Cat. No.	Description	Pk Qty
PS/96110	1ml pipette, single wrap	1000
PS/96118	5ml pipette, single wrap	200
PS/96128	10ml pipette, single wrap	200
PS/96134	25ml pipette, single wrap	200
PS/96310	1ml pipette, inner pack 25	1000
PS/96318	5ml pipette, inner pack 25	500
PS/96228	10ml pipette, inner pack 25	500

Pipette accessories

Rack, serological pipette, Acrylic

Storage rack for 1 to 100ml glass or plastic serological pipettes, 5 angled compartments with 2 holes on the top designed to house pipette controllers. Made from strong, clear Acrylic.

Cat. No.	Description
PP/12040	Rack, pipette, for 1 to 100ml pipettes, clear acrylic