

Pipette stand, vertical

Blue chemically resistant polypropylene stand. Very stable with 230mm diameter base plate. Holds up to 94 pipettes vertically.

Azlon

Cat. No.	Code	Description
PS/12034	SWP033	Pipette stand, vertical, holds 94




PS/12034

Pipette tray, draw or bench, horizontal, PVC

Can accommodate 7 pipettes with overall stem diameter not exceeding 20mm longitudinally or 16 pipettes with stem diameter not exceeding 10mm laterally. Very useful work station holder not only for pipettes.

Kartell

Cat. No.	Code	Dimensions (wxdxh), mm
PS/46208	996	280 x 215 x 40



PS/46208

Stand, pipette, polythene covered wire

Horizontal stand for up to 10 pipettes.

Cat. No.	Dimensions (wxdxh), mm
ST/93002	125 x 156 x 118mm




ST/93002

Pipette jar, PP

Rugged container for pipette storage. No more chipped ends as pipettes are cushioned by polypropylene. The jar, but not the PE cap can be autoclaved at 120°C for 20 minutes. Tight fitting cover protects pipettes during transportation and from dust etc. Welded flanged base gives excellent stability. Can be easily cut down to suit shorter length pipettes. Can be canted on side for drainage with pipettes in-situ.

Kartell

Cat. No.	Code	Diameter, mm	Height, inner/outer, mm
PJ/46202	592	82	410/426



PJ/46202