

Pipettors, single and multi-channel

Brand Transferpette® S pipettors feature a large central pipetting button to separate the ejection function and are suitable for both left and right single-handed operation.

Calibration is easy and does not require any special tools. The 4 position volume display is always clearly visible and an external flag indicates changes from factory settings. The pipettors can be completely autoclaved at 121°C (2bar) according to DIN EN 285

Supplied with DE-M marking, performance certificate, shelf/rack mount and silicone grease.

Transferpette S, fixed volume



PP/70002

Cat. No.	Code	Capacity
PP/70000	F-10	10µl
PP/70002	F-20*	20µl
PP/70004	F-25	25µl
PP/70006	F-50	50µl
PP/70008	F-100	100µl
PP/70010	F-200*	200µl
PP/70012	F-500	500µl
PP/70014	F-1000	1000µl

For use only with 20 - 200µl tips

Transferpette S, adjustable volume

Cat. No.	Code	Range
PP/70040	D-2.5	0.1 - 2.5µl
PP/70042	D-10	0.5 - 10µl
PP/70044	D-20	2 - 20µl
PP/70046	D-50	5 - 50µl
PP/70048	D-100	10 - 100µl
PP/70050	D-200	20 - 200µl
PP/70052	D-1000	100 - 1000µl
PP/70054	D-5000	500 - 5000µl
PP/70056	D-10000	1000 - 10000µl

Please note:

Model D-2.5, when using Brand nano-cap™ tips or tips from other manufacturers, calibration and possibly adjustments is necessary. Model D-5000 for use only with 5ml tip from Brand and 5ml tip from Thermo Fisher Scientific Inc. (Finnpipette®) D-10000 for use with 10ml tip from Brand and appropriate tips from Gilson® and Eppendorf®.

Pipette sets, Transferpette® S, adjustable volume in sets of 3 or 5

Cat. No.	Code	Description
PP/70080	MICRO	Pipette set of 3, D-2.5, D-10 and D-100 with 3 TipBoxes and 3 shelf/rack mounts
PP/70082	MIDI	Pipette set of 3, D-20, D-200 and D-1000 with 3 TipBoxes and 3 shelf/rack mounts
PP/70084	MACRO	Pipette set of 3, D-1000, D-5000 and D-10000 with 3 TipBoxes and 3 shelf/rack mounts
PP/70086	STANDARD	Pipette set of 3, D-10, D-100 and D-1000 with 3 TipBoxes and 3 shelf/rack mounts
PP/70088	PACKAGE 1	Pipette set of 5, D-2.5, D-10, D-100, D-200 and D-1000 with 3 TipBoxes and benchtop rack
PP/70090	PACKAGE 2	Pipette set of 5, D-10, D-100, D-200, D-1000 and D-5000 with 4 TipBoxes and benchtop rack

The Brand Transferpette® S-8 and -12 extends the range of single channel models to provide multi-channel pipetting and offers all the same features provided in the single channel range.

Supplied with DE-M marking, performance certificate, 2 x TipBox filled with tips, shelf/rack mount, reagent reservoir, mounting tool, silicone grease and 1 set of silicone/FKM sealing rings.

Transferpette S, multichannel

Cat. No.	Code	Number of channels	Range
PP/70100	M8-10	8	0.5 - 10µl
PP/70102	M8-50	8	5 - 50µl
PP/70104	M8-100	8	10 - 100µl
PP/70108	M8-300	8	30 - 300µl
PP/70120	M12-10	12	0.5 - 10µl
PP/70122	M12-50	12	5 - 50µl
PP/70124	M12-100	12	10 - 100µl
PP/70128	M12-300	12	30 - 300µl

Accessories:

Cat. No.	Code	Description
PS/70000	7048 05	Rack, bench-top to hold 6 single or multichannel pipettors
PS/70002	7048 10	Rack/shelf mount for one single or multichannel pipettor
PS/70004	7053 85	Stand, bench-top for single pipettor up to 1000µl capacity
PS/70006	7053 86	Stand, bench-top for single pipettor 500 to 5000µl capacity
PS/70008	7034 40	Stand, bench-top for single multichannel pipettor



PP/70100



PP/70120