



PT/37200

## Pipettor, multi-dispenser

Dispetta dispenser, robust and well-designed, hand operated, for serial dispensing using disposable tips. Five volume settings combined with seven dispenser tip sizes result in 35 different volumes from 10  $\mu$ l to 5 ml. The dispenser tip is filled similar to a syringe using a well-positioned lever, after preselecting the dispensing volume the pipet button is pressed for each dispensing step.

- Ergonomic design and light weight for convenient handling
- Smooth and fatigue-free operation
- Precise and convenient volume setting
- Clearly visible volume adjustment, optimized readability
- Table with volume scale fixed directly onto the instrument
- Precise gear mechanism for accurate and reproducible pipetting
- Safe, rapid and easy exchange of tips
- Maintenance free
- Robust design to ensure long life of unit

### Volume settings for tips

Capacity	1 $\mu$ l	2 $\mu$ l	3 $\mu$ l	4 $\mu$ l	5 $\mu$ l
0.5ml	10 $\mu$ l	20 $\mu$ l	30 $\mu$ l	40 $\mu$ l	50 $\mu$ l
1.25ml	25 $\mu$ l	50 $\mu$ l	75 $\mu$ l	100 $\mu$ l	125 $\mu$ l
2.5ml	50 $\mu$ l	100 $\mu$ l	150 $\mu$ l	200 $\mu$ l	250 $\mu$ l
5.0ml	0.25ml	0.5ml	0.75ml	1.0ml	1.25ml
12.5ml	0.5ml	1.0ml	1.5ml	2.0ml	2.5ml
25ml	1.0ml	2.0ml	3.0ml	4.0ml	5.0ml



PT/37220



PT/37222



PT/37224

Cat. No.	Description
PP/37200	Pipettor, Dispetta for serial dispensing, 5 $\mu$ l to 5ml
<b>Tips</b>	
PT/37220	Dispetta tip, non-sterile, 0.5ml, pack of 100
PT/37222	Dispetta tip, non-sterile, 1.25ml, pack of 100
PT/37224	Dispetta tip, non-sterile, 2.5ml, pack of 100
PT/37226	Dispetta tip, non-sterile, 5.0ml, pack of 100
PT/37228	Dispetta tip, non-sterile, 12.5ml, pack of 100
PT/37230	Dispetta tip, non-sterile, 25ml, pack of 50, requires adapter
PT/37232	Dispetta tip, non-sterile, 50ml, pack of 25, requires adapter
PT/37234	Adapter for use with 25 and 50ml non-sterile tips
PT/37250	Dispetta tip, sterile, 0.5ml, pack of 100
PT/37252	Dispetta tip, sterile, 1.25ml, pack of 100
PT/37254	Dispetta tip, sterile, 2.5ml, pack of 100
PT/37256	Dispetta tip, sterile, 5.0ml, pack of 100
PT/37258	Dispetta tip, sterile, 12.5ml, pack of 100
PT/37260	Dispetta tip, sterile, 25ml, pack of 25, requires adapter
PT/37262	Dispetta tip, sterile, 50ml, pack of 25, requires adapter
PT/37264	Adapter, sterile, for use with 25 and 50ml sterile tips



PT/37226



PT/37228



PT/37232



PT/37230



PT/37234