



PT/160



PT/37000



PT/37002



PT/37004



PT/37006



PT/37008



PT/37010



PT/37014

Pipette tips

Pipette tips, PP, non-sterile



Precision moulded with cleanly finished ends allowing good visibility of tip contents

- Manufactured from metal free virgin polypropylene
- Autoclavable tips are provided non-sterile and clean

Cat. No.	Code	Description	Volume μ l	Inner Pk	Pk Qty
PT/16002	BCT10	Clear Gilson Micro	0.5 – 10	1000	2000
PT/16004	BCT20	Clear Eppendorf Micro	0.5 – 10	1000	2000
PT/16006	BCT25	Yellow Gilson	5 – 200	1000	2000
PT/16008	BCT30	Yellow Universal/Eppendorf	2 – 200	1000	2000
PT/16016	BCT70	Blue Universal/Eppendorf	50 – 1000	1000	1000
PT/16106	BCB70	Universal Blue	50 - 1000	1000	10000
PT/16024	BCT100	Blue Gilson Macro	1000 – 5000	250	1000

Pipette tips, universal, non-sterile

Ratiolab® pipette tips are ideal for all commercially available pipettes. Their uniform high quality and precision make today's pipettes fully effective. The tips are suitable for use with manual or electronic, single or multichannel pipettes.

- Made of high quality polypropylene
- Tight fitting precision cone
- Autoclavable
- Guaranteed dust-free due to special ionization treatment
- Available pre-sterilized in Multiracks or non-sterile in bags
- Excellent transparency due to high quality surfaces
- Minimal residual wetting due to controlled quality of raw materials
- Excellent volume accuracy due to precise tip orifice
- Suitable for all commonly used pipette types (unless otherwise stated)

Cat. No.	Code	Description	Length mm	Pk Qty
PT/37000	2500190	Universal Pipette Tips crystal, 0,1-10 μ l, for Gilson pipetman/Neo	38.3	1000
PT/37002	2500170	Universal Pipette Tips crystal, 0.5-20 μ l	45.8	1000
PT/37004	2100600	Universal Pipette Tips yellow, 200 μ l	49.5	1000
PT/37006	2500160	Universal Pipette Tips yellow, 1-250 μ l	51.2	1000
PT/37008	2100610	Universal Pipette Tips blue, 100-1200 μ l	69.0	1000
PT/37010	2400650	Universal Pipette Tips Makro II, 1000-5000 μ l for Biohit Proline	133.3	300
PT/37012	2410650	Universal Pipette Tips Makro II, 1000-5000 μ l in rack for Biohit Proline	133.3	1 x 50
PT/37014	2400642	Universal Pipette Tips 1-10ml for Biohit mLine, Gilson Pipetman and Eppendorf Research Plus	153.6	100



PT/37012