

### Adapters, reduction

Cat. No.	Socket size	Cone size
DA01	10/19	14/23
DA12	14/23	19/26
DA13	14/23	24/29
DA14	14/23	29/32
DA23	19/26	24/29
DA24	19/26	29/32
DA25	19/26	34/35
DA34	24/29	29/32
DA35	24/29	34/35
DA36	24/29	40/38
DA37	24/29	45/40
DA38	24/29	50/42
DA39	24/29	55/44
DA45	29/32	34/35
DA47	29/32	45/40
DA56	34/35	40/38
DA57	34/35	45/40
DA58	34/35	50/42



DA14

### Adapters, expansion

Cat. No.	Socket size	Cone size
XA10	14/23	10/19
XA21	19/26	14/23
XA31	24/29	14/23
XA32	24/29	19/26
XA41	29/32	14/23
XA42	29/32	19/26
XA43	29/32	24/29
XA52	34/35	19/26
XA53	34/35	24/29
XA54	34/35	29/32



XA42

### Adapters, multiple, two necks

Parallel.

Cat. No.	Socket size	Cone size
MA1/11	14/23	14/23
MA1/2	19/26	19/26
MA1/3	19/26	24/29
MA1/33	24/29	24/29



MA1/3

### Adapters, multiple, three necks

Two parallel and one at 45°.

Cat. No.	Socket size	Cone size
MA4/23	19/26	24/29



MA4/23